This quick start guide is intended to provide a quick overview of your pipette's key features and to offer basic instructions for getting started. For detailed information, please refer to the operating instructions that can be found at www.integra-biosciences.com in different languages.

Intended Use

This is a general-purpose laboratory instrument for use in research only. VIAFLO/VOYAGER Electronic Pipettes are used for aspirating and dispensing aqueous liquids in the volume range of 0.5–5000 µl using GripTip™ pipette tips only, see www.integra-biosciences.com.

Battery Charge


Charge the battery completely before first use.


The supply voltage must meet the requirements of the mains adapter: 100–240 VAC, 50/60 Hz.

Low Battery

Blinking red


Recharging Battery

Blinking bars

~ 2.5 h

Fully Charged Battery

Three solid bars

Getting Started

1. **Battery Charge Indicator**
2. **Back Button**, to navigate backward
3. **Touch Wheel**, spin to scroll the cursor
4. **OK Button**, to make a selection
5. **Arrow Buttons**, for selections and tip spacing (VOYAGER only)
6. **Purge Button**, to empty tips
7. **Run Key**, to start operations
8. **Tip Ejector**

Turn on:

Press and release the Run Key.

Turn off:

Press and hold the Back Button for 3 seconds.

Language Selection

From the Main Menu, use the Touch Wheel to scroll to the Toolbox and press OK. Use the “Language” option to select the desired language.

Attaching GripTips

Press the pipette into the GripTips of the appropriate size until you hear and feel a click indicating that a seal has been achieved. When loading GripTips using a multichannel pipette, slowly rock the pipette from one side to the other side to ensure the proper seal is achieved.

Selecting a Pipetting Mode


With your electronic pipette, you have the choice of building protocols by selecting from a menu of predefined programming modes, or you can create custom multi-stepped protocols.

Pipetting mode	Description
Pipet	Liquid transfers when aspirate and dispense volumes are equal.
Repeat Dispense	Dispense multiple aliquots of the same volume.
Sample Dilute	Aspirate two liquids divided by an air gap, followed by dispense.
Pipet/Mix	Multiple mixing by aspiration and dispensing of defined volume.
Manual Pipet	Control the aspiration and dispensing up to the set volume.
Reverse Pipet	Liquid transfers of viscous or high vapor pressure liquids.
Variable Dispense	Dispense multiple aliquots of different volumes.
Multi Aspirate	Aspirate multiple aliquots of different volumes.
Sample Dilute/Mix	Aspirate two liquids divided by an air gap, followed by dispense and mix.
Serial Dilution	Aspirate a transfer volume followed by dispense and mix sequences.
Custom	Allows to create and store of up to 40 multi-stepped pipetting protocols.

From the Main Menu, use the Touch Wheel to scroll to your desired function-based programming mode. Press OK to enter the mode and to start defining parameters.


Setting/Changing Variables

Scroll to Edit and press OK.

A list of editable steps is displayed.

Select a step and press OK. Use the Touch Wheel to set the value and press OK. Press > to save your settings.

Pipette's Toolbox

The Toolbox provides a variety of options, as setting personal preferences like the Touch Wheel spin sensitivity, calibration, computer connectivity and device information.

Running a Program

Insert the GripTips into the liquid to be transferred. Press and release the Run Key to aspirate the volume selected in the first step (marked with a triangle on the Run screen).


Your pipette will prompt you to press the Run Key at each successive step.

Two Step Blowout

Note: During the last dispense of a program, you can perform a two step blowout to prevent liquid from being aspirated back into the GripTips:


- Press and hold the Run Key to dispense without blowin.
- Remove the GripTips from the target vessel.
- Release the Run Key to start blowin.

Tip Spacing (VOYAGER only)

In each pipetting mode, the available tip spacing positions are displayed at the bottom of the screen. The current position is highlighted.

At any program step press ◀ or ▶ to move the tip position to the next spacing.

To modify the tip spacing positions go to the “Tip Spacing” menu, accessible from the main menu of the pipette or within every pipetting mode.

Select the positions to be changed (First, Middle, Last) and press OK. Use ◀ and ▶ to move the tips to the desired position, click OK and then SAVE.

Equipment Disposal

Electronic pipettes must not be disposed of with unsorted municipal waste. Dispose of electronic pipettes in accordance with the laws and regulations in your area governing disposal of devices.

Manufacturer

INTEGRA Biosciences AG
CH-7205 Zizers, Switzerland
T +41 81 286 95 30
F +41 81 286 95 33

info@integra-biosciences.com
www.integra-biosciences.com

INTEGRA Biosciences Corp.
Hudson, NH 03051, USA
T +1 603 578 5800
F +1 603 577 5529

English (EN): Safety information

This is a general-purpose laboratory instrument for use in research only. Any use of this instrument in a medical or IVD setting is under the sole responsibility of the user.

- 1) Use only an INTEGRA battery and charging device.
- 2) Do not expose the battery to heat over 60 °C or mechanical stress.
- 3) Replace battery if the charging intervals are unusually short or if the charging takes much longer than 4 hours.
- 4) Do not use or charge the instrument in an atmosphere with danger of explosion.

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Declaration of conformity

INTEGRA Biosciences AG – 7205 Zizers, Switzerland

declares on its own responsibility that the devices

Description	Models
VIAFLO Pipettes	4011, 4012, 4013, 4014, 4015, 4016, 4621, 4622, 4623, 4624, 4626, 4631, 4632, 4633, 4634, 4636, 4641, 4642, 4646
VOYAGER Pipettes	4721, 4722, 4723, 4724, 4726, 4731, 4732, 4736, 4743, 4744, 4763, 4764

comply with:

EU Directives (DoW: Date of Withdrawal)	Before DoW	DoW	After DoW
Low Voltage Equipment	2006/95/EC	20.04.2016	2014/35/EU
Electromagnetic Compatibility	2004/108/EC	20.04.2016	2014/30/EU
Restriction of Hazardous Substances	2011/65/EU		
Waste Electrical and Electronic Equipment	2012/19/EU		
Battery Directive	2006/66/EC		

EU Regulations

Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH)	1907/2006
Capacity Labelling of Portable Secondary Batteries	1103/2010

Standards for EU

Safety requirements for electrical equipment for measurement, control and laboratory use - General requirements.	EN 61010-1: 2010
Electrical equipment for measurement, control and laboratory use - EMC requirements.	EN 61326-1: 2013

Standards for Canada and USA

Safety requirements for electrical equipment for measurement, control and laboratory use - General requirements.	CAN/CSA-C22.2 No. 61010-1
Safety requirements for electrical equipment for measurement, control and laboratory use - General requirements.	UL 61010-1
Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.	Part 15 of the FCC Rules Class A

International Standards

Piston-operated volumetric apparatus - Part 2: Piston pipettes	ISO 8655-2
--	-------------------

Zizers, March 12, 2018

 Elmar Morscher
 CEO

 Thomas Neher
 Quality Manager