



This quick start guide is intended to provide a quick overview of your ASSIST PLUS's key features and to offer basic instructions for getting started. For detailed information, please refer to the latest operating instructions (OI) that can be found at [www.integra-biosciences.com](http://www.integra-biosciences.com) in different languages.

### Intended use

This is a general-purpose laboratory instrument. Any use of this instrument in a medical or IVD setting is the sole responsibility of the user. ASSIST PLUS is used to perform pipetting operations automatically with a VIAFLO multichannel pipette, a VOYAGER pipette or a D-ONE single channel pipetting module attached.

### ⚠ Safety information

Regardless of the listed safety notes, all locally applicable regulations must be observed.

- 1) This product may only be operated in a secure, protected network with trustworthy clients.
- 2) The pipette may only be used by properly trained personnel in a manner specified by INTEGRA Biosciences.
- 3) Observe the hazard warnings on the device. Keep the tower moving area free of any labware and keep hands away.
- 4) Do not use or charge the instrument near flammable material or in an atmosphere with danger of explosion. Also, do not pipette highly flammable liquids such as acetone or ether.
- 5) Servicing work and repairs may only be performed by INTEGRA Biosciences or an authorized after-sales service member.

### Getting started



*Set up the instrument on a perfectly horizontal surface according to the IQ/OQ document (PN 128953). Supply voltage: 100–240 VAC, 47–63 Hz.*


Insert the adapter cable in the socket on the left side of the ASSIST PLUS and connect the adapter via a 3-core cable to the electricity mains.



- a. Gripper, with lever to release the pipette
- b. Tower, with pipette arm
- c. Interfaces and main switch
- d. Tip rack deck, for GRIPTIPS
- e. Touch panel
- f. Deck, with positions for labware
- g. Instrument base
- h. Waste bin, for tip ejection
- i. Tip sensor, for tip monitoring
- j. Reflector, for tip sensor

#### Switch on/off:

Press the main switch.



- k. Arrow keys to move the pipette back and fore (Y-axis), right and left (X-axis) and up and down (Z-axis)
- l. Start/pause key
- m. Communication LED and error LED
- n. On/standby key

# INTEGRA

Switch on the device. Press the blinking start/pause key (I) to home the ASSIST PLUS.



*Remove hands from ASSIST PLUS during homing and operation.*

## On the pipette: Installing the communication module

The VIAFLO/VOYAGER pipette or the D-ONE pipetting module can be connected to the ASSIST PLUS via the communication cable or the communication module, which must be installed:

1. Press and hold the Back button to turn off the device.



2. Loosen screw to detach pipette backing.
3. Disconnect battery power plug.



4. Insert the communication module into socket above battery. Position carefully and press to secure firmly.
5. Reconnect battery power plug and remount pipette back.

## On the pipette: Removing the tip ejector lid



For automatic tip ejection, remove the **tip ejector lid** with your thumb nail and store it in a safe place.

## On the pipette: Enabling ASSIST PLUS programs


The menu Toolbox – ASSIST PLUS provides a variety of options, like pairing ASSIST PLUS with a pipette, device information, tip monitoring, tip type, park position and more.



To display the ASSIST PLUS program in the main menu of the pipette, navigate to the Toolbox. Select Preferences - Main Menu and select ASSIST PLUS (green ✓).



## On the pipette: Pairing ASSIST PLUS

The first time a specific pipette is connected to ASSIST PLUS, both instruments need to be paired.



- Scroll to the Toolbox and press OK.
- From the ASSIST PLUS menu select ComModule Pairing.
- Switch OFF and ON the ASSIST PLUS device and wait until the message “Pairing successful” is displayed.


## Connection between ASSIST PLUS and pipette

When selecting ASSIST PLUS in the Main Menu, the pipette establishes a connection with the ASSIST PLUS. The communication symbol  is displayed next to the battery indicator and on the touch panel of the ASSIST PLUS the communication LED  (m) lights up in blue.

# INTEGRA

## Adapting the pipette holder

Adapt the pipette holder to the used INTEGRA multichannel pipette or D-ONE:



- VIAFLO 8- or 16-channel: fold both silver levers (1, 2) down.
- VIAFLO 12-channel: pull the lower silver lever (2) down.
- VOYAGER or D-ONE: fold both silver levers (1, 2) up.

## Attaching an electronic pipette

Rotate the lower housing of the pipette by 90 degrees as displayed in the picture below.



- To install a pipette, the tip ejector opening must face the pipette holder. Hold the pipette at an angle and insert it into the pipette holder. Then lift up until the black gripper snaps in place.
- To release the pipette, lift the black gripper by pushing down the lever.

## Setting up pipetting programs

You can easily create multi-step programs with the VIALAB software on a PC or select a predefined program on the pipette and adapt the parameters to your needs.

Program category	Description
VIALAB Programs	Contains all programs created with VIALAB software (the only option for D-ONE).
Serial Dilution	Aspirate a transfer volume followed by a mix.
Repeat Dispense	Dispense multiple aliquots of the same volume.
Variable Dispense	Dispense multiple aliquots of different volumes.
Multi Aspirate	Aspirate multiple aliquots of different volumes.
Plate Copy	Transfer samples between plates of the same well number.
Reformat	Reformat 12, 24, 48, 96 and 384 well plates.
Custom Programs	Create multi-stepped user-defined pipetting programs.

## Transfer VIALAB pipetting programs


Under Toolbox - Communications choose one of the following methods to transfer your VIALAB program to your pipette:

- 1) **USB:** Connect the pipette to the PC using the USB cable (#4226).
  - 2) **Via ASSIST PLUS:** Connect the pipette to ASSIST PLUS, then connect ASSIST PLUS and the PC via the USB cable (#4226) to transfer.
  - 3) **ComModule:** Connect the pipette and the PC via the Communication module (Bluetooth).
- Alternatively choose **PC control (A+)** to run your program directly from the PC:

## Running a program

When the pipette is connected to ASSIST PLUS, select from the ASSIST PLUS menu the program category and the program you would like to run. Follow the instructions on display. Place all required labware (reservoir, plates, etc.) on the deck in the correct orientation and position. D-ONE requires its own tip deck. No PC is required.

# INTEGRA



- Insert the tip box in the correct orientation and select the first row/column with tips. Press the RUN button of the pipette.
- Press the blinking start/pause key (I) on the ASSIST PLUS device. It switches to solid white; the ASSIST PLUS homes and the program is performed automatically.

A program can be paused by pressing the start/pause key (I). Either press the start/pause key again to continue the program or press Abort on the pipette to abort the program.

You can also run a program directly from your PC. Under Toolbox - Communications choose PC control (A+) and connect the pipette to ASSIST PLUS and ASSIST PLUS to the PC via USB cable. Once the connection is established, follow the instructions displayed in the VIALAB software.

## Maintenance



*Always switch off the device and disconnected from the electricity supply when carrying out maintenance work.*

- Clean the ASSIST PLUS housing with a lint-free cloth lightly soaked with mild soap solution in distilled water or with a 70 % solution of isopropanol or ethanol.
- To decontaminate the ASSIST PLUS, take out the waste bin (h) to the front, remove the waste bin frame, the magnetically attached sensor reflector (j) and the decks. Wipe the clean instrument surface with a lint-free cloth lightly soaked with a suitable disinfectant and wipe dry directly.
- Decontaminate the sensor reflector by spraying and wiping it with 70% Ethanol and after a short treatment time (1-2 minutes) wiping it dry.
- The aluminum decks and the waste bin with frame can be autoclaved at 121 °C for 20 min.
- The device may be decontaminated with H<sub>2</sub>O<sub>2</sub> gas (maximal concentration 35 %) for 60 minutes.

## Arm fixation for transport "Park Position"



- Remove the pipette, labware and deck.
- With the power on, press and hold the three buttons: fore, back and On.
- Once the arm lowers to the lowest central position and locks, turn the main switch OFF. The unit is now ready for transport.

## Equipment Disposal



ASSIST PLUS must not be disposed of with unsorted municipal waste. Dispose of ASSIST PLUS in accordance with the regulations in your area governing disposal of devices.



In certain regions and countries, e.g. in EU member states, the distributor is obliged to take back this product free of charge at the end of life. Please contact your local distributor for more details.

## Manufacturer and customer service

<b>INTEGRA Biosciences AG</b> Tardisstrasse 201 7205 Zizers, Switzerland	<b>INTEGRA Biosciences Corp.</b> 22 Friars Drive Hudson, NH 03051, USA	Find your local representative: <a href="http://www.integra-biosciences.com">www.integra-biosciences.com</a> <a href="mailto:info@integra-biosciences.com">info@integra-biosciences.com</a>
<b>Direct sales country</b>	<b>Telephone</b>	<b>E-mail</b>
Australia	+617 3497 5800	<a href="mailto:info-au@integra-biosciences.com">info-au@integra-biosciences.com</a>
Benelux	+31 630 609 866	<a href="mailto:info-benelux@integra-biosciences.com">info-benelux@integra-biosciences.com</a>
China	+86 21 5844 7203	<a href="mailto:info-cn@integra-biosciences.com">info-cn@integra-biosciences.com</a>
Denmark, Sweden, Norway	+45 3173 5373	<a href="mailto:info-nordic@integra-biosciences.com">info-nordic@integra-biosciences.com</a>
France	+33 1 34 30 76 76	<a href="mailto:info-fr@integra-biosciences.com">info-fr@integra-biosciences.com</a>
Germany, Austria	+49 6409 81999 15	<a href="mailto:info-de@integra-biosciences.com">info-de@integra-biosciences.com</a>
Ireland	+44 1635 797000	<a href="mailto:info-ie@integra-biosciences.com">info-ie@integra-biosciences.com</a>
Japan	+81 3 5962 4936	<a href="mailto:info-jp@integra-biosciences.com">info-jp@integra-biosciences.com</a>
Switzerland	+41 81 286 95 30	<a href="mailto:info-ch@integra-biosciences.com">info-ch@integra-biosciences.com</a>
United Kingdom	+44 1635 797000	<a href="mailto:info-uk@integra-biosciences.com">info-uk@integra-biosciences.com</a>
United States, Canada	+1 603 578 5800	<a href="mailto:info-us@integra-biosciences.com">info-us@integra-biosciences.com</a>



## Declaration of conformity

**INTEGRA Biosciences AG – 7205 Zizers, Switzerland**

declares on its own responsibility that the device

Description	Model
<b>ASSIST PLUS</b>	<b>4505</b>
<b>Mains adapter</b>	<b>128909</b>

complies with:

### EU Provisions

2014/35/EU, 2014/30/EU, 2014/53/EU, 2012/19/EC, 2011/65/EC, 1907/2006, 2019/1782

### UK Provisions

S.I. 2016/1101, S.I. 2016/1091, S.I. 2017/1206, S.I. 2013/3113, S.I. 2012/3032

Zizers, May 4, 2026

Urs Hartmann  
CEO

Daniela Gross  
Head of Corporate Quality

For detailed CE declaration with regulations of other countries, please refer to the operating instructions.