



This quick start guide is intended to provide a quick overview of your D-ONE single channel pipetting module's key features. For detailed information, please refer to the latest operating instructions of D-ONE and ASSIST PLUS that can be found at www.integra-biosciences.com in different languages.

Intended Use

This is a general-purpose laboratory instrument for use in research only. D-ONE is used on the ASSIST PLUS pipetting robot for aspirating and dispensing aqueous liquids in the volume range of 0.5–1250 µl using GripTip™ pipette tips only, see www.integra-biosciences.com. D-ONE provides pressure-based liquid level detection and is compatible with two different GripTip sizes.


⚠ Safety information

Regardless of the listed safety notes, all locally applicable regulations must be observed.

- 1) This product may only be operated in a secure, protected network with trustworthy clients.
- 2) The pipette may only be used by properly trained personnel in a manner specified by INTEGRA Biosciences.
- 3) Use only an INTEGRA battery and charging device.
- 4) Do not expose the battery to heat over 60 °C or mechanical stress.
- 5) Replace battery if the charging intervals are unusually short or if the charging takes much longer than 4 hours.
- 6) Do not use or charge the instrument near flammable material or in an atmosphere with danger of explosion.
- 7) Servicing work and repairs may only be performed by INTEGRA Biosciences or an authorized after-sales service member.



Battery Charge




Charge the Li-ion battery completely before first use with an INTEGRA stand, cable or on ASSIST PLUS. 



The supply voltage must meet the requirements of the mains adapter: 100–240 VAC, 50/60 Hz.



~ 2.5 h



Getting Started



1. **Battery charge indicator** on display
2. **Back button**, to navigate backward
3. **Touch wheel**, spin to scroll the cursor
4. **OK button**, to make a selection
5. **Arrow buttons**, for selections
6. **Purge button**, to empty tips
7. **RUN button**, to start operations
8. **Tip ejector**

Turn on:

Press and release **RUN**.

Turn off:

Press and hold the Back button for 3 seconds.

Language Selection

From the Main Menu, use the Touch Wheel to scroll to the Toolbox and press **OK**. Use the “Language” option to select the desired language.

Pipette’s Toolbox

The Toolbox provides a variety of options, as setting personal preferences like the Touch wheel spin sensitivity, computer connectivity and device information.

Creating a Pipetting Program

Use the VIALAB software to create multi-step programs on a PC and transfer them to the D-ONE, see VIALAB operating instructions.

Starting a Pipetting Program

Use the communication module or the communication cable to connect D-ONE with ASSIST PLUS.

The following GripTips for benchtop pipetting systems can be attached:

D-ONE 0.5–300 µl (# 4531): 12.5 µl (pink) and 300 µl (green).

D-ONE 5–1250 µl (# 4532): 125 µl (yellow) and 1250 µl (blue).

Please visit the GripTip Selector Guide to find the correct GripTips.



From the Main Menu, select ASSIST PLUS. Press **OK** to enter the list of VIALAB programs.

Select a program and press RUN.

Maintenance

Clean the outer surface with a lint-free cloth lightly soaked with mild soap solution in distilled water or with a 70 % dilution of Isopropyl or Ethanol. Never use Acetone or other solvents.

D-ONE cannot be disassembled or autoclaved. Decontaminate D-ONE in accordance with good laboratory practice. The device may be decontaminated with H₂O₂ gas (maximal concentration 35 %) for 60 minutes.

If you intend to ship D-ONE to be periodically calibrated, the original packaging can be stored and reused for this purpose.

Equipment Disposal



Electronic pipettes must not be disposed of with unsorted municipal waste.

Dispose of electronic pipettes in accordance with the regulations in your area governing disposal of devices.

In certain regions and countries, e.g. in all EU member states, the distributor is obliged to take back this product free of charge at the end of life. Please contact your local distributor for more details.

Manufacturer

INTEGRA Biosciences AG Tardisstrasse 201 CH-7205 Zizers, Switzerland T +41 81 286 95 30	info@integra-biosciences.com www.integra-biosciences.com	INTEGRA Biosciences Corp. 22 Friars Drive Hudson, NH 03051, USA T +1 603 578 5800
INTEGRA Biosciences SAS 8 avenue du Fief 95310 Saint Ouen l'Aumône, FR T +33 1 34 30 76 76	INTEGRA Biosciences Ltd 2 Rivermead Business Park Thatcham, Berks, RG19 4EP, UK T +44 1635 797 00	INTEGRA Biosciences KK Higashikanda 1-5-6, Chiyoda-ku Tokyo, 101-0031, JP T +813 5962 4936
INTEGRA Biosciences Deutschland GmbH An der Amtmannsmühle 1 35444 Biebertal, DE T +49 6409 81 999 15	INTEGRA Biosciences (Shanghai) Co., Ltd. Room 1110, No. 515 Huanke Road Shanghai 201315, CN T +86 21 5844 7203	



Declaration of conformity

INTEGRA Biosciences AG – 7205 Zizers, Switzerland

declares on its own responsibility that the devices

Description	Models
VIAFLO Pipettes	4011, 4012, 4013, 4014, 4015, 4016, 4621, 4622, 4623, 4624, 4626, 4631, 4632, 4633, 4634, 4636, 4641, 4642, 4646
VOYAGER Pipettes	4721, 4722, 4723, 4724, 4726, 4731, 4732, 4736, 4743, 4744, 4763, 4764
D-ONE Module	4531, 4532
Accessories	3214, 3215, 3216, 3217, 3218, 4200, 4205, 4210, 4211, 4215, 4221

comply with:

EU Directives and Regulations

2014/35/EU, 2014/30/EU, 2012/19/EC, 2011/65/EC, 2006/66/EC, 1907/2006, 2019/1782, 1103/2010

GBR Regulations

S.I. 2016/1101, S.I. 2016/1091, S.I. 2008/2164, S.I. 2013/3113, S.I. 2012/3032

Zizers, May 2, 2022

Urs Hartmann
CEO

Daniel Bächli
Head of Corporate Quality

For detailed CE declaration with regulations of other countries, please refer to the D-ONE operating instructions.