

# Minimizing dead volume during qPCR set-up with the divided reservoir on the ASSIST PLUS pipetting robot

## Introduction

Overcoming high dead volume to save precious and expensive reagents is an important goal for laboratories. The assumption that only pipetting from single tubes guarantees the lowest possible dead volume holds up well. However, even with automated single channel pipettes, this is not an option when high throughput is required.

INTEGRA has developed novel divided reservoirs to reduce reagent wastage when using the VOYAGER adjustable tip spacing pipette.

Reagent transfer can be fully automated with high accuracy and low dead volume by mounting the VOYAGER onto the ASSIST PLUS pipetting robot.

Here we demonstrate the cost- and time-saving benefits when preparing master mixes with expensive SYBR® Green dye, TaqMan® probes and DNA polymerase reagents for multiplex qPCR using a divided reservoir.

### Key benefits:

- Save expensive reagents like SYBR Green dye, TaqMan probes and DNA polymerase, and minimize the overall master mix dead volume down to 10 µl by using divided reagent reservoirs during automated multiplex qPCR set-up.
- The divided reagent reservoir combines the SureFlo™ anti-sealing array with a unique surface treatment that spreads liquid evenly, offering even lower dead volume.
- Hands-free time can be reduced by using a multichannel electronic pipette with a repeat dispense function, allowing multiple dispenses to be performed consecutively.
- The ASSIST PLUS pipetting robot guarantees accurate and fast reagent and sample transfer, helping to ensure the reproducibility of multiplex qPCR results.
- Up to four reagents can be processed simultaneously by using an ASSIST PLUS equipped with the VOYAGER adjustable tip spacing pipette and a dual reservoir adapter.
- INTEGRA's dual reservoir adapter and PCR cooling block can be chilled prior qPCR set-up, improving on-deck reagent stability for precious DNA polymerase.

## Overview: How to reduce dead volume during automated qPCR set-up


In this application note, we demonstrate how to carry out fast and accurate multiplex qPCR master mix transfers with the lowest possible dead volume using the ASSIST PLUS pipetting robot equipped with a VOYAGER adjustable tip spacing pipette and divided reagent reservoir.

When performing multiplex qPCR assays, SYBR Green dye, TaqMan probes targeting different genes of interest, a DNA polymerase and buffer are combined into the master mix, which is then transferred to a PCR plate. The method of transferring the master mix plays a crucial role in decreasing reagent waste and costs during automated qPCR set-up.



The divided reservoir is split into two compartments, with maximum volumes of 5 ml and 10 ml. The 5 ml compartment offers the lowest available dead volume, and can be accessed with the VOYAGER pipette by adjusting the tip spacing to the minimal distance. The SureFlo bottom allows aspiration from the very bottom of the reservoir to reduce the dead volume.

### Experimental set-up



**Figure 1:** Deck set-up for qPCR master mix transfer from the divided reagent reservoir. **Position A:** source – dual reservoir adapter with INTEGRA's divided reagent reservoir and qPCR master mix in the 5 ml compartment A1 (blue). **Position B:** target - Sapphire 96 well PCR plate on INTEGRA's cooling block (pink). **Position C:** empty.

The 8 channel 12.5 µl VOYAGER adjustable tip spacing pipette with 12.5 µl sterile, filter, low retention GRIPTIPS is mounted on the ASSIST PLUS pipetting robot.


On deck position A (**Figure 1**), INTEGRA's divided reagent reservoir with SureFlo anti-sealing array is set up onto the dual reservoir adapter with qPCR master mix in the 5 ml compartment, A1. An INTEGRA PCR cooling block together with a Sapphire 96 well PCR plate is placed in landscape orientation on deck position B (**Figure 1**).

### Step by step procedure:

Select and run the VIALAB program 'Master\_mix\_transfer\_from\_reservoir'. The 8 channel 12.5 µl VOYAGER adjustable tip spacing pipette aspirates 10 µl of the master mix from A1 of the divided reagent reservoir (**Figure 2**), and dispenses it into each well of the 96 well PCR plate. A pre- and post-dispense of 0.5 µl guarantees precision while pipetting expensive reagent.

### Tips:

- Pipetting height is crucial for accessing the last microliters of liquid. Set the Safety Bottom Offset to 0 mm, and pipette from the very bottom of the reservoir.
- Setting a slow aspiration speed for the last aspiration step ensures that the correct volume of liquid will aspirate when the liquid level is low.
- When using low retention GRIPTIPS, pre-wetting and pre-/post-dispense steps improve accuracy and precision while pipetting small volumes.
- If focusing on higher throughput over low dead volume, the VOYAGER adjustable tip spacing pipette can be exchanged for a higher volume model for repeat dispensing.



**Figure 2:** The dual reservoir adapter on the ASSIST PLUS pipetting robot accommodates two divided reservoirs.

## Results

To highlight the time and cost savings of INTEGRA's divided reservoir when used with the VOYAGER adjustable tip spacing pipette on the ASSIST PLUS pipetting robot, we investigated how quickly filling a 96 well PCR plate with 10  $\mu$ l of multiplex qPCR master mix could be achieved.


As shown in **Figure 3**, we compared the residual dead volume and hands-on times when transferring multiplex qPCR master mix in single dispenses (**Figure 3a**) or multiple dispenses (**Figure 3b**) with:

1. a single channel pipette from a microcentrifuge tube
2. a VOYAGER multichannel pipette from a standard or divided reservoir
3. a VOYAGER pipette mounted on an ASSIST PLUS from a divided reservoir


As there are no comparable reservoirs on the market available, we chose a 25 ml reservoir and used a similar set-up as for the divided reagent reservoir during master mix transfer.

Using a single channel manual pipette, the 96 well plate was filled with 10  $\mu$ l master mix in each well in 6 minutes, by individual transfers from a microcentrifuge tube filled with 970  $\mu$ l of master mix. The excess volume needed was reduced to 1 % of the total volume of the master mix (1 reaction volume). The processing time was decreased by using a multichannel pipette and a reservoir, especially when using the repeat dispense function of an electronic pipette. However, the use of the standard reservoir led also to an increase in the dead volume up to 260  $\mu$ l (26 reaction volumes), which equates to 23 % more master mix needed per 96 well plate compared to the microcentrifuge tube. When using the divided reagent reservoir and the VOYAGER instead, this issue was avoided completely, resulting in only a 1 % dead volume (1 reaction volume). This matched the advantage of using microcentrifuge tubes, while keeping the speed benefits of a multichannel pipette.

a)




b)



**Figure 3:** Illustration of dead volume and operating time after filling twenty 96 well qPCR plates (a) without and (b) with repeat dispensing.

Automating the VOYAGER on the ASSIST PLUS using INTEGRA's divided reagent reservoir led to additional hands-free time, and throughput was further increased by using the repeat dispensing function. Automation provided 2 hours of hands-free time while filling twenty 96 well qPCR plates with master mix, and successfully kept the dead volume down to 1 % (20 reaction volumes). Compared to the standard reservoir, the divided reservoir saved 480 reaction volumes (**Figure 4**). This emphasizes the potential cost savings of the reservoirs, particularly for performing multiplex qPCR with SYBR Green, large amounts of expensive TaqMan probes, and precious DNA polymerase.



**Figure 4:** More than 20 % qPCR reaction volumes saved with INTEGRA's divided reservoir.

## Remarks

- **VIALAB software:** The VIALAB programs can be easily adapted to your specific labware and protocols, for instance when partial plates are needed.
- **Partial plates:** Programs can be adapted at any time to a different number of samples, giving laboratories total flexibility to meet current and future demands.

## Conclusion

- It is possible to decrease the dead volume of expensive multiplex qPCR master mix down to 1 % when using INTEGRA's divided reagent reservoir and VOYAGER adjustable tip spacing pipette.
- Pipetting reagents from the divided reservoir with a VOYAGER pipette is much faster, and achieves the same low dead volume as pipetting from a tube with a single channel pipette.
- The ASSIST PLUS pipetting robot equipped with the VOYAGER is the optimal set-up when pipetting small volumes of expensive reagents like SYBR Green, TaqMan probes and DNA polymerase, enabling slow aspiration speeds and no interference from shaky hands.
- Thanks to INTEGRA's dual reservoir adapter and divided reagent reservoirs, up to four different reagents can be processed under cooled conditions at the same time.

## Materials

Manufacturer	Part Number	Description	Link
INTEGRA Biosciences	4505	ASSIST PLUS base unit	<a href="https://www.integra-biosciences.com/global/en/pipetting-robots/d-one-for-assist-plus">https://www.integra-biosciences.com/global/en/pipetting-robots/d-one-for-assist-plus</a>
INTEGRA Biosciences	4721	VOYAGER 8 channel 12.5 µl electronic pipette	<a href="https://www.integra-biosciences.com/global/en/electronic-pipettes/voyager">https://www.integra-biosciences.com/global/en/electronic-pipettes/voyager</a>
INTEGRA Biosciences	4547	Dual Reservoir Adapter	<a href="https://www.integra-biosciences.com/global/en/pipetting-robots/assist-plus">https://www.integra-biosciences.com/global/en/pipetting-robots/assist-plus</a>
INTEGRA Biosciences	6250	PCR Cooling Block 96	<a href="https://www.integra-biosciences.com/global/en/pipetting-robots/assist-plus">https://www.integra-biosciences.com/global/en/pipetting-robots/assist-plus</a>
INTEGRA Biosciences	6505	12,5 µl LONG, Sterile, Filter Low Retention GRIPTIPS	<a href="https://www.integra-biosciences.com/global/en/pipette-tips/griptip-selector-guide">https://www.integra-biosciences.com/global/en/pipette-tips/griptip-selector-guide</a>
INTEGRA Biosciences	4356	25 ml Divided Reservoir, Polypropylene	<a href="https://www.integra-biosciences.com/global/en/reagent-reservoirs/multichannel-reagent-reservoirs">https://www.integra-biosciences.com/global/en/reagent-reservoirs/multichannel-reagent-reservoirs</a>
Greiner Bio-One International	652270	96 Well PCR Sapphire Plate	<a href="https://shop.gbo.com/en/switzerland/products/bio-science/molecular-biology/pcr-microplates/652270.html?sword_list%5B0%5D=652270&amp;no_cache=1&amp;ga=2.40164706.1588166689.1664192269.799770341.1664192267">https://shop.gbo.com/en/switzerland/products/bio-science/molecular-biology/pcr-microplates/652270.html?sword_list%5B0%5D=652270&amp;no_cache=1&amp;ga=2.40164706.1588166689.1664192269.799770341.1664192267</a>

**INTEGRA Biosciences AG**  
7205 Zizers, Switzerland  
T +41 81 286 95 30  
F +41 81 286 95 33  
info@integra-biosciences.com

**INTEGRA Biosciences Corp.**  
Hudson, NH 03051, USA  
T +1 603 578 5800  
F +1 603 577 5529  
info-us@integra-biosciences.com

**INTEGRA Biosciences Deutschland GmbH**  
35444 Biebertal, Deutschland  
T +49 6409 81 999 15  
F +49 6409 81 999 68  
info-de@integra-biosciences.com

**INTEGRA Biosciences SAS**  
95062 Cergy-Pontoise Cedex 1, France  
T +33 (0)1 34 30 76 76  
F +33 (0)1 34 30 76 79  
info-fr@integra-biosciences.com

**INTEGRA Biosciences Ltd.**  
Thatcham, Berks RG19 4EP, UK  
T: +44 1635 797000  
F: +44 1635 797001  
info-uk@integra-biosciences.com

**INTEGRA Biosciences Nordic ApS**  
Vallensbækvej 22A 3TV  
Brøndby 2605, Denmark  
T + 45 3173 5373  
info-nordic@integra-biosciences.com

**INTEGRA Biosciences (Shanghai) Co., Ltd.**  
中国上海自由贸易试验区环科路515号1110室  
邮编: 201315  
电话: +86 21 5844 7203  
info-cn@integra-biosciences.com

インテグラ・バイオサイエンス株式会社  
〒101-0031  
東京都千代田区東神田 1-5-6  
東神田MK第五ビル 3階  
T 03-5962-4936  
F 03-5962-4937  
info-jp@integra-biosciences.com