

High-throughput, Three-Dimensional Assay Development with Corning® 1536-well Spheroid Microplates and INTEGRA Biosciences' VIAFLO 384 and VIAFILL Liquid Handlers

CORNING

Application Note

Hilary Sherman¹, Marc Hamel², Francesca Chang², and Ann Rossi², Ph.D.

¹Corning Incorporated, Life Sciences, Kennebunk, ME USA

²INTEGRA Biosciences, Hudson, NH USA

Introduction

It is becoming more apparent that 3D models often better recapitulate the *in vivo* cell environment than traditional 2D methodologies. Accordingly, more spheroid and other 3D models are being utilized for drug discovery. To meet the demands of most high-volume drug discovery programs, a high throughput system is essential. The Corning 1536-well spheroid microplate can generate 1536 uniform, single spheroids that can be assayed via imaging, fluorescence, or luminescence directly in the microplate. With a maximum volume of 14 μL per well, accurate liquid handling is essential to the generation of quality data. We have established an automated method for seeding, dosing, and assaying spheroids in a Corning 1536-well spheroid microplate format with the assistance of INTEGRA Biosciences' liquid handling instruments. The VIAFILL is designed for rapid bulk liquid dispensing at volumes as low as 0.5 μL and can be fully automated with a plate stacker accessory. The VIAFLO 384 is a handheld electronic pipetting system which enables transfer into 24-, 96-, 384-, and 1536-well microplates in a fast, compact, and easy-to-use manner. Here we have demonstrated robust quality data for assay development and high throughput screens generated via this method.

Materials and Methods

Spheroid Formation

DU 145 (ATCC® HTB-81) and PANC-1 (ATCC CRL-1469™) cells were routinely cultured in Dulbecco's Modification of Eagle's Medium (DMEM; Corning Cat. No. 10-013-CM) containing 10% fetal bovine serum (FBS; Corning Cat. No. 35-010-CV). Cells were seeded into 1536-well spheroid microplates (Corning Cat. No. 4527) at 1×10^3 cells per well in a volume of 5 μL per well using the VIAFILL with

the 16-channel dispensing cassette (INTEGRA Biosciences Cat. No. 5742). Microplates were incubated overnight at 37°C in a humidified CO₂ incubator to allow for spheroid formation.

Homogenous Cell-based Assay

After DU 145 spheroid formation, 1 μL of 6 μM staurosporine (MilliporeSigma Cat. No. 569396) or buffer control was dispensed into each well via the VIAFLO 384 with a 384-channel pipetting head (0.5-12.5 μL) and 12.5 μL short GripTips (INTEGRA Biosciences Cat. No. 6475). The next day, 6 μL of CellTiter-GLO® 3D (Promega Cat. No. G9683) was added to each well and allowed to incubate at room temperature for 1 hour. Luminescent signal was then read via a TECAN Infinite® M200 microplate reader.

High-content Imaging Assay

After spheroid formation, both DU 145 and PANC-1 spheroids were exposed to various concentrations of cisplatin (Sigma-Aldrich/MilliporeSigma Cat. No. 1134357) by adding 1 μL via the VIAFLO 384 with 12.5 μL short GripTips. The plates were incubated overnight. To assess cytotoxicity, 6 μL per well of 40 $\mu\text{g}/\text{mL}$ Hoechst 34580 (Thermo Fisher Cat. No. H21486) and 16 $\mu\text{g}/\text{mL}$ propidium iodide (PI; AnaSpec Cat. No. 83215) in phosphate buffered saline (PBS; Corning Cat. No. 21-030-CM) were added to each well for 1 hour. Microplates were then imaged on the Thermo Fisher CellInsight™ CX7 high-content screening platform.

Results

Assay Robustness

The unique design of the Corning 1536-well spheroid microplate allows for the generation of a single, uniform spheroid in each well (Figure 1). To ensure that the spheroids are of uniform size




Figure 1. Uniform, single spheroids formed in each well across the entire spheroid microplate. Representative image of single uniform DU 145 spheroids formed in the Corning 1536-well spheroid microplate with one well digitally zoomed in. Images were taken with Thermo Fisher CellInsight CX7 high-content screening platform using a 4X objective.



Figure 2. Microplate dispenser for uniform cell seeding. INTEGRA VIAFILL for rapid, bulk liquid dispensing into 96-, 384-, or 1536-well microplates.

across the entire microplate, accurate liquid handling of the cell suspension during seeding is critical. For rapid cell seeding, the INTEGRA VIAFILL was utilized (Figure 2). This bulk liquid dispenser can accurately seed a 1536-well microplate in less than 1 minute per plate and can be paired with a plate stacker for increased throughput. For compound and reagent addition, the INTEGRA VIAFLO 384 was utilized. With its small footprint, the handheld electronic 24-, 96-, and 384-channel pipetting system can increase throughput and productivity in filling microplates while using disposable tips (Figure 3). The uniformity of signal, signal-to-background ratio (S/B), and calculated Z' values were assessed to demonstrate the utility of this experimental design for 3D cell screening. Uniformity of signal was demonstrated by treating sections of the microplate with a single concentration of buffer or staurosporine to induce cell death and assaying spheroid viability via homogenous luminescence cell viability assay. A heat map of luminescence signal shows good signal uniformity within




Figure 3. Electronic pipetting system for accurate compound delivery. INTEGRA VIAFLO 384 for increased productivity and efficiency to accurately work with 24-, 96-, 384-, and 1536-well microplates.

treatment conditions regardless of location in the microplate (Figure 4). The mean S/B ratio from this assay was greater than 20 for 3 independent studies (Figure 5). Additionally, Z' values were consistently above 0.7, indicating the robustness of the viability assay in the 1536-well spheroid microplate (Figure 5).

Proof of Concept Assay

DU 145 and PANC-1 spheroids exposed to varying concentrations of cisplatin were stained and imaged with the Thermo Fisher CellInsight™ CX7 high-content screening platform. Representative images were captured directly without the need to transfer the spheroids to another microplate, demonstrating utility of culture and subsequent imaging of spheroids in the spheroid microplate (Figure 6). Image analysis of PI staining reveals typical concentration-response curves for cisplatin cytotoxicity with IC₅₀ values of 0.03375 and 0.03917 mM for DU 145 and PANC-1 spheroids, respectively (Figure 7).






Figure 5. Large, robust luminescence assay window. (A) Mean \pm SD luminescence of DU 145 spheroids exposed to buffer or 1 μ M staurosporine, resulting in S/B of 20.2. $p < 0.0001$ with unpaired t-test. (B) Excellent Z' values above 0.7 in 3 independent studies. N = 2304 wells in 3 independent studies.




Figure 6. Confocal imaging of spheroids within Corning 1536-well spheroid microplate. Representative z-stack images of DU 145 (top) and PANC-1 (bottom) spheroids exposed to 0 mM, 0.025 mM, and 83 mM cisplatin (left, middle, right, respectively) for 24 hours. Spheroids were stained with Hoechst (blue) and PI (red) to assess cell viability. Spheroids were imaged directly in the spheroid microplate with the Thermo Fisher CellInsight™ CX7 high-content screening platform using a 10X objective. Scale bar = 100 μ m.




Figure 7. Concentration-dependent cisplatin cytotoxicity. Mean \pm SD PI fluorescence of DU 145 and PANC-1 spheroids following 24 hours exposure to varying concentrations of cisplatin. N = 96 wells.

Conclusions

To screen large compound libraries against 3D cell models the proper tools are necessary to achieve quality results. The Corning® 1536-well spheroid microplate allows for the formation of uniform, single spheroids in each well that are ideally suited for high throughput 3D drug screens. The black sidewalls and clear bottom design of the Corning spheroid microplate make it ideally suited for imaging, luminescent, and fluorescent cell-based assays. Due to the small well volume of 1536-well microplates, accurate and consistent dispensing are essential. The INTEGRA VIAFILL and VIAFLO 384 instruments utilized for liquid handling for the 1536-well spheroid microplate result in assays with high signal-to-background ratios with excellent and consistent Z' values.

For more specific information on claims, visit the Certificates page at www.corning.com/lifesciences.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit www.corning.com/lifesciences or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

CORNING

Corning Incorporated
Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

ASIA/PACIFIC

Australia/New Zealand
t 61 427286832

China
t 86 21 3338 4338
f 86 21 3338 4300

India
t 91 124 4604000
f 91 124 4604099

Japan

t 81 3-3586 1996
f 81 3-3586 1291

Korea

t 82 2-796-9500
f 82 2-796-9300

Singapore

t 65 6572-9740
f 65 6735-2913

Taiwan

t 886 2-2716-0338
f 886 2-2516-7500

EUROPE

CSEurope@corning.com

France

t 0800 916 882
f 0800 918 636

Germany

t 0800 101 1153
f 0800 101 2427

The Netherlands

t 020 655 79 28
f 020 659 76 73

United Kingdom

t 0800 376 8660
f 0800 279 1117

All Other European Countries

t +31 (0) 206 59 60 51
f +31 (0) 206 59 76 73

LATIN AMERICA

grupoLA@corning.com

Brasil

t 55 (11) 3089-7400

Mexico

t (52-81) 8158-8400