

VIAFLO 96 | VIAFLO 384

24, 96 And 384 Channel Handheld Electronic Pipette
Easy, Productive, And Affordable

VIAFLO 96 | VIAFLO 384


Handheld microplate processing

Have you reached the limits of productivity with traditional handheld pipettes? Would you like to increase your pipetting throughput and reproducibility without investing in large and expensive automation?

VIAFLO 96 and VIAFLO 384 enable transfer of 24, 96, and 384 samples in a single step, thus increasing pipetting throughput and reproducibility significantly. At a fraction of the cost, they often surpass throughput capabilities of expensive and complicated robots.

VIAFLO 96

VIAFLO 96 is a handheld 24 or 96 channel pipette designed to fit on any lab bench. The unique operating concept makes it as easy to use as a traditional handheld pipette.



VIAFLO 384

VIAFLO 384 is an enhanced version of the successful VIAFLO 96 pipette. It offers the same intuitive operation and compact design, but allows working with 24, 96 and 384 channel heads for an even higher productivity.

Intuitive use

VIAFLO 96 and VIAFLO 384 offer the ease of use of traditional handheld pipettes. No position programming or special training is required. To move the pipetting head to the desired position, simply push in that direction. The servo assisted steering mechanism supports the movement for an effortless workflow and precise tip positioning.



A choice of ten predefined pipetting modes provides users with convenient solutions for the most common liquid handling tasks.

Standard pipetting routines, such as plate to plate transfers or multiple dispenses from a reservoir to different plates, can be executed with minimal set up time.


For more elaborate pipetting routines, up to 40 personalized custom programs can be created and saved.

All settings of the pipette are easily modified over the touch wheel controlled graphical user interface.

Accelerate multiwell pipetting


VIAFLO 96 and 384 offer unmatched speed when transferring 24, 96, or 384 samples at once. It is ideal for fast and precise sample transfers either from plate to plate, reservoir to plate or within plates.

VIAFLO 96 is **4 times faster** than an 8 channel electronic pipette for filling five 96 well plates. Additionally, the risk for pipetting errors is minimized, which provides better accuracy and precision.



Reservoir to plate transfers


- Compound addition
- Medium exchange
- Cell seeding
- Addition of buffers and reagents
- 24, 96, 384, 1536 well



The „Repeat Dispense“ mode allows for a very efficient filling of multiple microplates with the same liquid.

Plate to plate transfers

- Plate duplication
- Supernatant transfers/collections
- Plate reformatting
- 24 to 96 well, 96 to 384 well, 384 to 1536 well



An exact aspiration height can be defined to precisely remove supernatant and ensures pipette tips do not impact the bottom of a plate.

Sample transfers within plates

- Serial dilutions
- Mixing and resuspension



Performing serial dilutions and mixing of entire plates is easily accomplished. The number of mixing cycles and mixing speed can be exactly defined.

Flexible volume range

Users can simply change the pipetting head of their VIAFLO 96 and VIAFLO 384 to guarantee optimal matching of the available volume range to the application performed. Changing a pipetting head of VIAFLO 96 and VIAFLO 384 is easy and takes less than one minute. Each channel of a pipetting head is independent, having its own precise piston and sealing to exclude any cross-channel influences for error-free pipetting.





Plate holder options


VIAFLO 96 and VIAFLO 384 feature two deck positions (A and B) to accommodate different plate holders. They can be removed and replaced by another option any time to adapt the instrument to various application requirements.



Standard plate holder

Plate holder with slide function

For 24 and 96 channel pipetting. The standard plate holder (6205) is used for simple 24 or 96 well transfers. The plate holder with slide function (6210) enables a grid offset to access 384 well plates with a 96 channel head.



Spring loaded plate holder with slide function

For 24, 96 and 384 channel pipetting. Available for position A (6215) and B (6220), these plate holders feature a spring load mechanism to secure microplates in place. Recommended for applications that require an exact plate positioning. Both holders also feature the slide function.

Plate holder inserts: PCR cooling blocks

These aluminum 96 and 384 well PCR cooling blocks keep reagents chilled during sample preparation in a PCR plate. The cooling blocks feature a SBS footprint, fit into any INTEGRA plate holders and can also be used for other applications or with other instruments.

Three position stage

For 96 and 384 channel pipetting. The stage provides three plate holder positions. Move the stage to the left and right to work with all three positions. It features the spring mechanism and slide function

on all three positions. The three position stage occupies both deck positions A and B and cannot be combined with another plate holder.

Visit our website to find more plate holder options in the VIAFLO 96/384 parts and numbers section.

Technical data

	VIAFLO 96	VIAFLO 384
Channels	24, 96	24, 96, 384
Pipetting technology	air displacement	
Compatible plate formats	24, 96, 384, and 1536 well microplates, shallow and deep well	
Dimensions (W x D x H)	42 cm x 30 cm x 54 cm (17" x 12" x 21")	
Weight	25.7 kg (56.7 lbs)	27.7 kg (61.1 lbs)
Operating conditions	Temperature: 5 °C – 40 °C, Humidity: < 85 % non-condensing	
Power requirements	100 – 240 VAC, 50/ 60 Hz	

Head specifications

Part no.	6101 / 6131	6106 / 6136	6102 / 6132	6103 / 6123	6104 / 6124
Vol. Range	0.5 – 12.5 µl	2 – 50 µl	5 – 125 µl	10 – 300 µl	50 – 1250 µl
Accuracy	± 1.5 %	± 1.5 %	± 1.5 %	± 1.2 %	± 1.8 %
Precision	≤ 1.0 %	≤ 1.0 %	≤ 1.0 %	≤ 0.8 %	≤ 0.8 %

Specifications presented for nominal volume. Please refer to our website for a complete list of specifications.



Small space requirements

With its compact design, this benchtop system requires little space on a lab bench and fits into most laminar flow hoods to allow work under sterile conditions. Up to three stage positions accommodate tip racks, microplates, reagent reservoirs or other SBS-size labware. Up to 8 microplates can be stacked to make maximum use of the stage positions.


GripTips for reliable pipetting results

Reliable tip sealing and firm tip attachment to the pipette are essential for reproducible and error-free pipetting.

VIAFLO 96 and VIAFLO 384 features the advanced GripTip system. During loading the GripTip pipette tips are locked onto the tip fitting of the pipetting head, which guarantees that they will never accidentally fall off, even when performing a tip touch on the sidewall of a vessel. In addition, the system provides a perfect sealing of every tip attached, preventing uneven pipetting or sample spillage.

The GripTips are comfortably loaded and ejected without any physical effort due to a power assisted mechanism activated by a push button. The loading procedure also guarantees that GripTips are not touched by the user and remain sterile.

GripTips are offered in non-sterile, sterile and sterile filtered configurations to serve all application requirements. Filtered GripTips reliably protect the pipetting head and samples against contaminations caused by aerosols.



GripTips are packaged in two different tip racks dependent on the pipette tip volume. Larger volume GripTips, 300 μ l and 1250 μ l, are packaged in racks of 96 tips. Smaller volume GripTips, 12.5 μ l and 125 μ l, are packaged in space efficient racks of 384 tips. These XYZ racks contain four times as many tips and are molded from an anti-static material to prevent static build up.

GripTips for use with VIAFLO 96 and VIAFLO 384 undergo a strict quality control to assure perfectly straight pipette tips, enabling accurate dispenses.

Ordering information

Visit the *parts and numbers* section on our website to view all available accessories and consumables.

1. Choose a base instrument

Description	Part no.
VIAFLO 96	6001
24 and 96 channel handheld pipette	
VIAFLO 384	6031
24, 96, and 384 channel handheld pipette	

2. Select one or more pipetting heads

Description	Volume	Part no.	Compatibility	
			VIAFLO 96	VIAFLO 384
24 channel pipetting heads	10 – 300 µl	6123	✓	✓
	50 – 1250 µl	6124	✓	✓
96 channel pipetting heads	0.5 – 12.5 µl	6101	✓	✓
	2 – 50 µl	6106	✓	✓
	5 – 125 µl	6102	✓	✓
	10 – 300 µl	6103	✓	✓
	50 – 1250 µl	6104	✓	✓
384 channel pipetting heads	0.5 – 12.5 µl	6131	-	✓
	2 – 50 µl	6136	-	✓
	5 – 125 µl	6132	-	✓

3. Add one plate holder for each deck position on A and B

For more plate holder options, please refer to our website.

Description	Part no.
Standard Plate Holder – For 24 and 96 channel applications	6205
Plate Holder with slide function – For plate reformatting 96 well ↔ 384 well	6210
Spring Loaded Plate Holder A with slide function – Recommended for 384 well applications and plate reformatting	6215
Spring Loaded Plate Holder B with slide function – Recommended for 384 well applications and plate reformatting	6220
Three Position Stage (occupies both deck positions A and B) – Suitable for 96 and 384 well applications	6230

GripTips for use with VIAFLO 96 and VIAFLO 384

For more GripTip options, including Low Retention, Wide Bore, and long/short versions, please refer to our website.

Tip Size	Description	Filter Sterile	Sterile	Non-sterile
		Part no.	Part no.	Part no.
12.5 µl	5 XYZ racks of 384 tips	6455	6454	6453
50 µl / 125 µl	5 XYZ racks of 384 tips	6465	6464	6463
300 µl	5 racks of 96 tips	6435	6434	6433
	GREENCHOICE Refills, 5 inserts of 96 tips	-	-	6432
1250 µl	5 racks of 96 tips	6445	6444	6443
	GREENCHOICE Refills, 5 inserts of 96 tips	-	-	6442

Automation Friendly Reagent Reservoirs

Description	Part no.
Reservoir Bases	150 ml, 8 Reservoir Bases, Non-Sterile, Automation Friendly (microplate format)
	300 ml, 8 Reservoir Bases, Non-Sterile, Automation Friendly (microplate format)
Reservoir Inserts	150 ml, 30 Individually Sealed, Sterile, incl. 1 base and 1 lid
	150 ml, 100 Bulk, Sterile, incl. 1 base and 1 lid
	300 ml, 30 Individually Sealed, Sterile, incl. 1 base and 1 lid
	300 ml, 100 Bulk, Sterile, incl. 1 base and 1 lid
Reservoir Lids	150 ml, 25 Lids, Sterile
	300 ml, 25 Lids, Sterile

Accessories

Description	Part no.
PCR 96 well Cooling Block	6250
PCR 384 well Cooling Block	6255

INTEGRA Biosciences AG
7205 Zizers, Switzerland
T +41 81 286 95 30
F +41 81 286 95 33
info@integra-biosciences.com

INTEGRA Biosciences Corp.
Hudson, NH 03051, USA
T +1 603 578 5800
F +1 603 577 5529
info-us@integra-biosciences.com

INTEGRA Biosciences Ltd.
Thatcham, Berks RG19 4EP, UK
T: +44 1635 797000
F: +44 1635 797001
info-uk@integra-biosciences.com



www.integra-biosciences.com

INTEGRA