This quick start guide is intended to provide a quick overview of your pipette's key features and to offer basic instructions for getting started. For detailed information, please refer to the latest operating instructions that can be found at www.integra-biosciences.com in different languages.

Intended Use

This is a general-purpose laboratory instrument for use in research only. VIAFLO/VOYAGER Electronic Pipettes are used for aspirating and dispensing aqueous liquids in the volume range of 0.5–5000 µl using GripTip™ pipette tips only, see www.integra-biosciences.com.

Battery Charge


Charge the battery completely before first use.


The supply voltage must meet the requirements of the mains adapter: 100–240 VAC, 50/60 Hz.

Low Battery


Recharging Battery

~ 2.5 h

Fully Charged Battery

Getting Started

1. **Battery Charge Indicator**
2. **Back Button**, to navigate backward
3. **Touch Wheel**, spin to scroll the cursor
4. **OK Button**, to make a selection
5. **Arrow Buttons**, for selections and tip spacing (VOYAGER only)
6. **Purge Button**, to empty tips
7. **Run Key**, to start operations
8. **Tip Ejector**

Turn on:

Press and release the Run Key.

Turn off:

Press and hold the Back Button for 3 seconds.

Language Selection

From the Main Menu, use the Touch Wheel to scroll to the Toolbox and press OK. Use the “Language” option to select the desired language.


Attaching GripTips

Press the pipette into the GripTips of the appropriate size until you hear and feel a click indicating that a seal has been achieved. When loading GripTips using a multichannel pipette, slowly rock the pipette from one side to the other side to ensure the proper seal is achieved.

Selecting a Pipetting Mode


With your electronic pipette, you have the choice of building protocols by selecting from a menu of predefined programming modes, or you can create custom multi-stepped protocols.

Pipetting mode	Description
Pipet	Liquid transfers when aspirate and dispense volumes are equal.
Repeat Dispense	Dispense multiple aliquots of the same volume.
Sample Dilute	Aspirate two liquids divided by an air gap, followed by dispense.
Pipet/Mix	Multiple mixing by aspiration and dispensing of defined volume.
Manual Pipet	Control the aspiration and dispensing up to the set volume.
Reverse Pipet	Liquid transfers of viscous or high vapor pressure liquids.
Variable Dispense	Dispense multiple aliquots of different volumes.
Multi Aspirate	Aspirate multiple aliquots of different volumes.
Sample Dilute/Mix	Aspirate two liquids divided by an air gap, followed by dispense and mix.
Serial Dilution	Aspirate a transfer volume followed by dispense and mix sequences.
Custom	Allows to create and store of up to 40 multi-stepped pipetting protocols.

From the Main Menu, use the Touch Wheel to scroll to your desired function-based programming mode. Press OK to enter the mode and to start defining parameters.


Setting/Changing Variables

Scroll to Edit and press OK.

A list of editable steps is displayed.


Select a step and press OK. Use the Touch Wheel to set the value and press OK. Press > to save your settings.

Pipette's Toolbox

The Toolbox provides a variety of options, as setting personal preferences like the Touch Wheel spin sensitivity, calibration, computer connectivity and device information.


Running a Program

Insert the GripTips into the liquid to be transferred. Press and release the Run Key to aspirate the volume selected in the first step (marked with a triangle on the Run screen).

Your pipette will prompt you to press the Run Key at each successive step.

Two Step Blowout

During aspiration, the piston of your electronic pipette moves up (1).

During dispensing, the piston returns to the initial position (0). During the last dispense of a program, it automatically moves further down (2) and blows the remaining liquid from the tip (**Blowout**↓).

When the piston returns to the initial position (0), a small amount of air is aspirated, provided the tip is no longer immersed in the liquid (**Blowin**↑).

Note: Manually delay the blowin by holding the Run Key pressed during the last dispense. Remove the tips from the vessel and release the Run Key to start blowin.

Tip Spacing (VOYAGER only)

In each pipetting mode, the available tip spacing positions are displayed at the bottom of the screen. The current position is highlighted.

At any program step press ◀ or ▶ to move the tip position to the next spacing.

To modify the tip spacing positions go to the “Tip Spacing” menu, accessible from the main menu of the pipette or within every pipetting mode.

Select the positions to be changed (First, Middle, Last) and press OK. Use ◀ and ▶ to move the tips to the desired position, click OK and then SAVE.

Equipment Disposal

Electronic pipettes must not be disposed of with unsorted municipal waste.

Dispose of electronic pipettes in accordance with the laws and regulations in your area governing disposal of devices.

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English (EN): Safety information

This is a general-purpose laboratory instrument for use in research only. Any use of this instrument in a medical or IVD setting is under the sole responsibility of the user.

- 1) Use only an INTEGRA battery and charging device.
- 2) Do not expose the battery to heat over 60 °C or mechanical stress.
- 3) Replace battery if the charging intervals are unusually short or if the charging takes much longer than 4 hours.
- 4) Do not use or charge the instrument in an atmosphere with danger of explosion.

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INTEGRA Biosciences AG – 7205 Zizers, Switzerland