

---

PRESS RELEASE

**Electronic Pipettor Offers Affordable and Flexible Sample Handling...**

**INTEGRA** has announced the introduction of the **Voyager™ adjustable spacing pipettor**

Available with six, eight or twelve channels the new Voyager™ electronic pipettors feature a revolutionary mechanism that alters tip spacing at the press of a button. Voyager incorporates a small motor that moves the center-to-center tip spacing quickly and smoothly at any point in a protocol, all with one hand. Users can preset up to three distinct tip spacings to accommodate almost any labware vessels. Voyager Tip spacing is adjustable from 4.5mm to 19.5mm.

All Voyager pipettors share the same easy-to-use Touch Wheel™ controller from the award winning Viaflo® pipettors, with intuitive programming and the industry's first full-color display. Also, all Voyager pipettors are designed to utilize Viaflo Griptips. As a result, individuals can pipette without ever having a tip loosen and leak or fall off. A total of eight Voyager models are now available in volume ranges from 0.5ul to 1250ul. This unique product now allows users to easily work between 24, 48, 96, 384 cell culture and assay plates, tube racks, sample storage gel electrophoresis boxes and many additional labware formats.

INTEGRA is a leading provider of high-quality laboratory tools for liquid handling, media preparation, sterilization and cell cultivation. The entire organization is committed to creating innovative solutions which fulfil the needs of our customers in research, diagnostics and quality control within the life science markets and medical industry.

INTEGRA, comprising INTEGRA Biosciences AG and the Viaflo Corporation, maintains a worldwide network of over 90 independent highly trained distribution partners. INTEGRA Biosciences AG is headquartered in Zizers, Switzerland; Viaflo Corporation is located in Hudson, NH, USA.

For further information on the productivity gains achievable with Voyager please contact INTEGRA Biosciences AG on telephone +41-81-286-9530 or by email [info@integra-biosciences.com](mailto:info@integra-biosciences.com). In the USA please contact Viaflo Corporation on telephone +1-603-578-5800 or by email [loneil@viaflo.com](mailto:loneil@viaflo.com).

Illustrative image: (image available on request)



**JANUARY 2010**

**ibspr109.doc**

For further information: Dr Bill Bradbury (tel. +44-208-546-0869 / email [info@primetek-solutions.com](mailto:info@primetek-solutions.com))