

Kanal	Art.-Nr.	Volumen-Bereich (µl)	Volumen-schritt (µl)	Test-Volumen (µl)	Genauigkeit (±%)	Präzision (≤%)
24	6123	10 - 300	0.5	30	2.50	1.50
				150	1.50	1.00
				300	1.20	0.80
	6124	50 - 1250	1	125	2.50	1.50
				625	2.00	1.00
				1250	1.80	0.80
96	6101	0.5 - 12.5	0.01	1.25	8.00	3.00
				6.25	2.00	1.50
				12.5	1.50	1.00
	6106	2 - 50	0.05	5	5.00	2.50
				25	2.00	1.50
				50	1.50	1.00
	6102	5 - 125	0.1	12.5	3.00	1.70
				62.5	2.00	1.50
				125	1.50	1.00
	6103	10 - 300	0.5	30	2.50	1.50
				150	1.50	1.00
				300	1.20	0.80
	6104	50 - 1250	1	125	2.50	1.50
				625	2.00	1.00
				1250	1.80	0.80
384	6131	0.5 - 12.5	0.01	1.25	8.00	3.00
				6.25	2.00	1.50
				12.5	1.50	1.00
	6136	2 - 50	0.05	5	5.00	2.50
				25	2.00	1.50
				50	1.50	1.00
6132	5 - 125	0.1	12.5	3.00	1.70	
			62.5	2.00	1.50	
			125	1.50	1.00	

1. Die Spezifikationen bilden die Eigenschaften aller Pipettierkanäle ab.
2. Spezifikationen gelten nur für Einzelpipettierungen im Modus "Pipettieren".
3. Präzision = Variationskoeffizient