

关于符合中国《电子信息产品污染控制管理办法》的声明 宣言 INTEGRA FIREBOY

China RoHS declaration for INTEGRA FIREBOY



INTEGRA Biosciences maintains a direct sales and support organization in the United States of America, France, Germany, Austria and Switzerland, as well as a network of over 100 trained distribution partners worldwide












Relevant environment information is posted on the company's web site www.integra-biosciences.com and the correct URL for the web site is informed to the users through operating instructions and markings on the instruments.

I. Environmental Protection Use Period (EPUP)

EPUP is defined as a period in which, when toxic and hazardous substances or elements contained in the electronic information products are in normal use, such substances or elements will not leak out or undergo abrupt change, and the use of the electronic information product by an electronic information product user will not result in serious pollution to the environment or result in serious personal injury or property damage with regard to the user.

EPUP for INTEGRA Biosciences FIREBOY products is declared in the following table

Instruments	Description	Part #	EPUP
	FIREBOY plus, FIREBOY eco	144000, 144010	

Accessories	Description	Part #	EPUP
	Mains adapters	156630-156634	
	Battery NiMH 3.6V/2000 mAh	144310	
	Foot switch and connecting cable	143200	
	Benchtop switch	171081	
	Long burner head for heating in combination with e.g. tripod	144306	
	Windshield for hoods	143205	
	Borosilicate glass splash protector	143210	
	Gas cartridge adapters	144050-144054	
	Gas safety tubing	140055	
	Adapters	144227, 144240	

For inquiries on the date of manufacture please contact your local distributor or INTEGRA Biosciences AG by mail info@integra-biosciences.com or use the information request form on www.integra-biosciences.com. State your request and make sure you provide the serial number of the product of interest.

II. Hazardous Substance Content Disclosure Table

The disclosure is a binary yes/no disclosure that indicates that one homogeneous material within the subassembly according to the EIP-A, B, and C definitions in SJ/T 11363-2006 is exceeding the maximum concentration values.

关于符合中国《电子信息产品污染控制管理办法》的声明 宣言 INTEGRA FIREBOY

Management Methods on Control of Pollution from Electronic Information Products Declaration for INTEGRA FIREBOY

产品中有毒有害物质的名称及含量

Name and content of hazardous substances in products

这个表被制备，根据 SJ / T11364-2014 的规定

This table is prepared in accordance with the provisions of SJ/T11364-2014

分总成 Subassembly	零件号 Part Number	有毒有害物质或元素 Hazardous Substances or Elements					
		铅 (Pb)	汞 (Hg)	镉 (Cd)	六价铬 (Cr6+)	多溴联苯 (PBB)	多溴二苯醚 (PBDE)
壳体上部 Case upper part	144000, 144010	O	O	O	O	O	O
操作按钮 Knobs		O	O	O	O	O	O
阀块 Valve		O	O	O	O	O	O
燃烧器单元 Burner unit		O	O	O	O	O	O
液晶显示器 (选修的) Display (optional)		O	O	O	O	O	O
电路板 Electronics		O	O	O	O	O	O
电池 (选修的) Battery (optional)		O	O	O	O	O	O
电源 Mains adapter	156630-156634	X	O	O	O	O	O
附件 Accessories	144310, 143200, 171081, 144306, 143205, 143210, 144050- 144054, 140055, 144227, 144240	O	O	O	O	O	O

O: 表示该有毒有害物质在该部件所有均质材料中的含量均在 SJ/T 11363-2006 标准规定的限量要求 以下。

O: Indicates that this hazardous substance contained in all homogeneous materials of this part is below the limit requirement in SJ/T 11363-2006

X: 表示该有毒有害物质至少在该部件的某一均质材料中的含量超出 SJ/T 11363-2006 标准规定的限 量要求。

X: Indicates that this hazardous substance contained in at least one of the homogeneous materials of this part is above the limit requirement in SJ/T 11363-2006.

对销售之日的所售产品,本表显示我公司供应链的电子产品信息产品可能包含这些物质。注意：在所售产 品中可能 会也可能不会含有所有所列的部件。

This table shows where these substances may be found in the supply chain of our electronic information products, as of the date of sale of the enclosed product. Note that some of the component types listed above may or may not be a part of the enclosed product.