

TITLE: **VIAFLO / VOYAGER Generations**
 AUTHOR: TBe, 3.4.18

Over the course of time INTEGRA’s pipettes received updates to their hardware to add new features. This resulted in different pipette Generations.

Firmwares are not compatible with instruments between different Generations. For example, Generation 2 firmware is not compatible with Generation 3 instruments, and vice versa.

	Hardware Generation	Firmware	Serial number
VIAFLO / VOYAGER	Gen 1	1.00 - 2.99	below 600000
	Gen 2	3.00 - 3.99	6000001 - 6999999
	Gen 3	4.00 and higher	7000001 and higher

For VIALINK compatibility with the different firmware versions, see *VIAFLO Software and firmware update history*. This can be accessed directly from VIALINK

