

# Maximize Your Handheld Pipetting Productivity!


## Affordable VIAFLO 96 Bundle

### Easy to Use

- Unique operating concept mimics a handheld pipette.
- Enables transfers of 96 and 384 samples in a single step, significantly increasing pipetting throughput and reproducibility.

### Small Footprint

- With a compact design, this benchtop system requires little space on a lab bench and fits into most laminar flow hoods, ensuring a sterile environment.



| Part | Part Description                                   | Includes | Bundle Price      |
|------|--|----------|-------------------|
| 6001 | VIAFLO 96 Base Unit                                | 1        | <h1>\$14,995</h1> |
| 6101 | 96 Channel Pipetting Head 0.5-12.5 µl              | Choose 1 |                   |
| 6106 | 96 Channel Pipetting Head 2-50 µl                  |          |                   |
| 6102 | 96 Channel Pipetting Head 5-125 µl                 |          |                   |
| 6103 | 96 Channel Pipetting Head 10-300 µl                |          |                   |
| 6104 | 96 Channel Pipetting Head 50-1250 µl               |          |                   |
| XXXX | Any Case of 6000 Series GripTips                   | 1        |                   |
| XXXX | Any Case of Automation Friendly Reagent Reservoirs | 1        |                   |
| 6205 | Standard Plate Holder                              | 1        |                   |
| 6210 | Plate Holder For 384 Well Plate Indexing           | 1        |                   |



The standard plate holder (6205) is used for simple 96-well transfers. The plate holder with slide function (6210) enables a grid offset to access 384-well plates with a 96-channel head.



6205



6210

To benefit from this offer please contact INTEGRA Biosciences by [email](#) or phone and mention the promotion code **VIAFLO 96 Bundle** or clearly write promotion code on all faxed purchase orders. This offer is valid until April 30, 2018.

**Offer valid in the United States only.**

INTEGRA Biosciences, US  
2 Wentworth Drive  
Hudson, NH 03051

Phone: 1-603-578-5800  
Fax: 1-603-577-5529  
E-Mail: [us@integra-biosciences.com](mailto:us@integra-biosciences.com)

[www.integra-biosciences.com](http://www.integra-biosciences.com)

**INTEGRA**