



America

CERTIFICATE


No. U8 10 01 42035 005

Holder of Certificate: **INTEGRA Biosciences AG**

Tardisstrasse 201
7205 Zizers
SWITZERLAND

Production Facility(ies): 42035

Certification Mark:



Product: **Laboratory Equipment (media preparation)**

Model(s): **MEDIACLAVE 10**

Parameters:

Rated voltage:	200-240V
Rated frequency:	50/60Hz
Rated power consumption:	3750W
Protection class:	I

Tested according to:

- UL 61010-1:2004
- IEC 61010-2-010:2003
- CAN/CSA-C22.2 No. 61010-1:2004
- CAN/CSA-C22.2 No. 61010-2-010:2004

The product was voluntarily tested according to the relevant safety requirements and mentioned properties. It can be marked with the certification mark shown above. The certification mark must not be altered in any way. See also notes overleaf.

Test report no.: 028-71338437-000

Date, 2010-01-12



Page 1 of 1