



CERTIFICATE

No. U8 10 01 42035 006

Holder of Certificate: **INTEGRA Biosciences AG**

Tardisstrasse 201
7205 Zizers
SWITZERLAND

Production Facility(ies): 42035

Certification Mark:



Product: **Laboratory Equipment (media preparation)**

Model(s): **MEDIACLAVE 30**

Parameters:	Rated voltage:	3x200-208V	3x400V
	Rated frequency:	50/60Hz	50/60Hz
	Rated power consumption:	9300W	9300W
	Protection class:	I	

Tested according to: UL 61010-1:2004
IEC 61010-2-010:2003
CAN/CSA-C22.2 No. 61010-1:2004
CAN/CSA-C22.2 No. 61010-2-010:2004

The product was voluntarily tested according to the relevant safety requirements and mentioned properties. It can be marked with the certification mark shown above. The certification mark must not be altered in any way. See also notes overleaf.

Test report no.: 028-71338438-000

Date, 2010-01-12

Page 1 of 1

