

Save On CLEAR ADVANTAGE™ Reagent Reservoirs

Maximize reagent savings with the lowest dead volume reservoirs in the world!

1 to 2 cases – 35% off
3 or more cases – 50% off



SureFlo™ Multichannel Reagent Reservoirs

- SureFlo™ array spreads reagent evenly across the reservoir and allows tips to sit on the bottom without aspirating air. This further reduces dead volume, saving you money on reagents and maximizing sample recovery!
- Disposable inserts, with clearly visible graduations fit into reusable bases to reduce plastic waste.
- Convenient pour back spouts allow easy return of excess fluid to a source container, minimizing reagent waste.




SureFlo™ anti-sealing array


Dead Volume	Volume		Description	Polystyrene Part No.	List Price
< 20 µl	10 ml For use with 8 and 16 channel pipettes		30 Individually Sealed, Sterile, SureFlo	4371	CHF 64
			200 Bulk, Sterile, SureFlo	4372	CHF 134
			10 Reservoir Bases	4306	CHF 31
< 50 µl	25 ml For use with 8, 12 and 16 channel pipettes		30 Individually Sealed, Sterile, SureFlo	4381	CHF 68
			200 Bulk, Sterile, SureFlo	4382	CHF 141
			10 Reservoir Bases	4304	CHF 31
< 50 µl	100 ml For use with 8, 12 and 16 channel pipettes		30 Individually Sealed, Sterile, SureFlo	4391	CHF 72
			200 Bulk, Sterile, SureFlo	4392	CHF 146
			10 Reservoir Bases	4305	CHF 31

Divided Reagent Reservoir

- 5 and 10 ml compartments minimize dead volumes when working with smaller reagent volumes using multichannel pipettes or fill to a max volume of 25ml with a single fluid.
- The SureFlo™ array allows tips to rest directly on the bottom of the flat-bottom reservoir without seal off, for an incredibly low dead volume of less than 10 µl.





SureFlo™ anti-sealing array

Dead Volume	Volume	Description	Polystyrene Part No.	List Price	Polypropylene Part No.	List Price
< 10 µl small side < 20 µl large side	 25 ml	30 Individually Sealed, Sterile, SureFlo	4351	CHF 100	4356	CHF 129
		200 Bulk, Sterile, SureFlo	4352	CHF 208	4357	CHF 267
		10 Reservoir Bases	4304	CHF 31		

Automation Friendly Reagent Reservoirs

- Clear reservoirs in a reusable base, with integrated volume graduations for more accurate measurements and less wasted reagent.
- Dual viewing windows allow tips to be safely positioned at the bottom of the reservoir.
- Flat-bottom reservoir features a 96/384 SureFlo™ anti-sealing pattern to prevent pipette tip seal off and liquid ‘popping’ up into tips, filters or the pipetting head.
- Latching lid with an integrated gasket prevents spillage and evaporation during transport or storage.
- Formatted to fit in all SBS holders.



Dead Volume	Volume	Description	Polystyrene Part No.	List Price	Polypropylene Part No.	List Price
< 3 ml	 150 ml	8 Reservoir Bases, Non-Sterile	6301	CHF 36	6301	CHF 36
		25 Lids, Sterile	6302	CHF 51	6302	CHF 51
		30 Individually Sealed, Sterile, SureFlo	6317	CHF 112	6337	CHF 170
		100 Bulk, Sterile, SureFlo	6318	CHF 163	6338	CHF 246
< 3 ml	 300 ml	8 Reservoir Bases, Non-Sterile	6305	CHF 41	6305	CHF 41
		25 Lids, Sterile	6306	CHF 51	6306	CHF 51
		30 Individually Sealed, Sterile, SureFlo	6327	CHF 123	6347	CHF 187
		100 Bulk, Sterile, SureFlo	6328	CHF 180	6348	CHF 272

To benefit from this offer please contact INTEGRA Biosciences by **email** or phone and mention the promotion code **SureFlo Reservoirs** or clearly write promotion code on all faxed purchase orders. This offer is valid until February 28 2020.

Offer valid in Switzerland only.

INTEGRA Biosciences AG
7205 Zizers, Switzerland

T +41 81 286 95 55
F +41 81 286 95 07
info-ch@integra-biosciences.com

www.integra-biosciences.com

INTEGRA