

# Save On CLEAR ADVANTAGE™ Reagent Reservoirs

Maximize reagent savings with the lowest  
dead volume reservoirs in the world!

1 to 2 cases – 35% off  
3 or more cases – 50% off






## SureFlo™ Multichannel Reagent Reservoirs

- SureFlo™ array spreads reagent evenly across the reservoir and allows tips to sit on the bottom without aspirating air. This further reduces dead volume, saving you money on reagents and maximizing sample recovery!
- Disposable inserts, with clearly visible graduations fit into reusable bases to reduce plastic waste.
- Convenient pour back spouts allow easy return of excess fluid to a source container, minimizing reagent waste.




SureFlo™ anti-sealing array

Dead Volume	Volume		Description	Polystyrene Part No.	List Price
< 20 µl	10 ml For use with 8 and 16 channel pipettes		30 Individually Sealed, Sterile, SureFlo	4371	\$66
			200 Bulk, Sterile, SureFlo	4372	\$138
			10 Reservoir Bases	4306	\$26
< 50 µl	25 ml For use with 8, 12 and 16 channel pipettes		30 Individually Sealed, Sterile, SureFlo	4381	\$70
			200 Bulk, Sterile, SureFlo	4382	\$145
			10 Reservoir Bases	4304	\$26
< 50 µl	100 ml For use with 8, 12 and 16 channel pipettes		30 Individually Sealed, Sterile, SureFlo	4391	\$74
			200 Bulk, Sterile, SureFlo	4392	\$150
			10 Reservoir Bases	4305	\$26

## Divided Reagent Reservoir

- 5 and 10 ml compartments minimize dead volumes when working with smaller reagent volumes using multichannel pipettes or fill to a max volume of 25ml with a single fluid.
- The SureFlo™ array allows tips to rest directly on the bottom of the flat-bottom reservoir without seal off, for an incredibly low dead volume of less than 10 µl.




SureFlo™ anti-sealing array

Dead Volume	Volume	Description	Polystyrene Part No.	List Price	Polypropylene Part No.	List Price
< 10 µl small side < 20 µl large side	25 ml	30 Individually Sealed, Sterile, SureFlo	4351	\$103	4356	\$132
		200 Bulk, Sterile, SureFlo	4352	\$214	4357	\$275
		10 Reservoir Bases	4304	\$26		

## Automation Friendly Reagent Reservoirs

- Clear reservoirs in a reusable base, with integrated volume graduations for more accurate measurements and less wasted reagent.
- Dual viewing windows allow tips to be safely positioned at the bottom of the reservoir.
- Flat-bottom reservoir features a 96/384 SureFlo™ anti-sealing pattern to prevent pipette tip seal off and liquid 'popping' up into tips, filters or the pipetting head.
- Latching lid with an integrated gasket prevents spillage and evaporation during transport or storage.
- Formatted to fit in all SBS holders.



Dead Volume	Volume	Description	Polystyrene Part No.	List Price	Polypropylene Part No.	List Price
< 3 ml	150 ml	8 Reservoir Bases, Non-Sterile	6301	\$33	6301	\$33
		25 Lids, Sterile	6302	\$46	6302	\$46
		30 Individually Sealed, Sterile, SureFlo	6317	\$102	6337	\$175
		100 Bulk, Sterile, SureFlo	6318	\$148	6338	\$254
< 3 ml	300 ml	8 Reservoir Bases, Non-Sterile	6305	\$37	6305	\$37
		25 Lids, Sterile	6306	\$46	6306	\$46
		30 Individually Sealed, Sterile, SureFlo	6327	\$112	6347	\$193
		100 Bulk, Sterile, SureFlo	6328	\$163	6348	\$280

To benefit from this offer please contact INTEGRA Biosciences by [email](#) or phone and mention the promotion code **SureFlo Reservoirs** or clearly write promotion code on all faxed purchase orders. This offer is valid until December 31, 2019.

**Offer valid in the United States.**

INTEGRA Biosciences, US  
22 Friars Drive  
Hudson, NH 03051

Phone: 1-603-578-5800  
Fax: 1-603-577-5529  
E-Mail: [us@integra-biosciences.com](mailto:us@integra-biosciences.com)

[www.integra-biosciences.com](http://www.integra-biosciences.com)

**INTEGRA**