

CERTIFICATE OF COMPLIANCE FILTER STERILE Part No. 153015

This is to certify that the Class VI Biosafe product described below has been tested and released based on the results of approved test methods. The quality management system under which it was manufactured is in compliance with ISO 9001 and ISO 13485.

Description:	Filter 0.45 µm sterile
Part Number:	153015
Lot Number:	451037
Date of sterilization:	2010-09
Shelf life of sterility:	3 years
Expiration date:	2013-09
Housing material:	PVC, outlet is blue
Connectors:	Male slip luer
Filter material:	PTFE
Pore size rating:	0.45 µm
Packaging:	Single blister, sterile
Sterilization method:	EO (Ethylene Oxide)

BSE/TSE Declaration:

The filter product is produced with raw materials manufactured in conformity with the European Community Directive 2003/32/EC.

Thomas Neher
Quality Manager



2010-09-28

This certificate applies only to products in undamaged condition, which have been shipped, stored and / or used in appropriate manner.