

Improve the safety of your laboratory and save up to 30%!

VACUSAFE and **VACUSIP** liquid waste vacuum systems are the perfect replacement to in-house and self-made vacuum systems.

Both systems combine a vacuum pump, collection bottle, overflow protection and filter in an all-in-one liquid aspiration system. These systems offer the most effective method for containing biohazardous liquids with optimal personal protection.



SAVE UP TO 30%!

Which system is right for me?



VACUSAFE



VACUSIP

| | VACUSAFE | VACUSIP |
|--|--|--|
| Bottle Size | Polypropylene, 4 L/safety coated glass, 2 L | Glass, 0.5 L |
| Aspiration Range | 15 ml/s | 2.3 ml/s |
| Volumes for Aspiration | 2 ml - 500 ml | Up to 10 ml |
| Typical Labware | T-flasks, falcon tubes, petri dishes, multiwell culture plates, etc. | Microcentrifuge tubes, microplates, small falcon tubes, small petri dishes, etc. |
| Level Detection (stops aspiration before overfilling) | Yes (optional) | No |

| Part No. | Item Description | List Price | 1 Unit 10% Discount | 2-3 Units 20% Discount | 4+ Units 30% Discount |
|----------|---|-------------|------------------------|---------------------------|--------------------------|
| 158300 | VACUSAFE 2 L bottle (glass) with bottle handle and closed lid Bottle adapter Level detection and lid with quick tubing couplings (PVDF) Tubing set with filter VACUBOY hand operator (155510) | AUD \$4,676 | AUD \$4,208 | AUD \$3,741 | AUD \$3,273 |
| 158310 | VACUSAFE Polypropylene bottle 4 L with lid Lid for Polypropylene bottle with quick couplings and level detection Tubing set with filter Mains adapter: 100-240VAC, 50/60Hz VACUBOY hand operator (155510) | AUD \$3,666 | AUD \$3,299 | AUD \$2,933 | AUD \$2,566 |
| 158320 | VACUSAFE Polypropylene bottle 4 L with lid Lid for Polypropylene bottle with barbed tube fittings Tubing set with filter Mains adapter: 100-240VAC, 50/60Hz VACUBOY hand operator (155510) | AUD \$2,760 | AUD \$2,484 | AUD \$2,208 | AUD \$1,932 |
| 159000 | VACUSIP Aspiration System with battery | AUD \$1,265 | AUD \$1,139 | AUD \$1,012 | AUD \$886 |
| 159010 | VACUSIP Aspiration System without battery | AUD \$1,059 | AUD \$953 | AUD \$847 | AUD \$741 |



To benefit from this offer, please contact INTEGRA Biosciences by email or phone and mention the promotional code **Vacuum2024**, or clearly write the promotion code on faxed purchase orders. Discounts cannot be combined with any other offers. This offer is valid from May 1, 2024 through July 31, 2024.

INTEGRA Biosciences Pty Ltd
Unit 55, 193-203 South Pine Road
Brendale QLD 4500
Australia

Phone: +61 7 3497 5800
E-Mail: orders-au@integra-biosciences.com

www.integra-biosciences.com



Offer valid in Australia only.