

Improve the safety of your laboratory and save up to 30%!

VACUSAFE and **VACUSIP** liquid waste vacuum systems are the perfect replacement to in-house and self-made vacuum systems.

Both systems combine a vacuum pump, collection bottle, overflow protection and filter in an all-in-one liquid aspiration system. These systems offer the most effective method for containing biohazardous liquids with optimal personal protection.



Which system is right for me?



VACUSAFE



VACUSIP

Bottle Size
Aspiration Range
Volumes for Aspiration
Typical Labware
Level Detection (stops aspiration before overflowing)

Polypropylene, 4 L/safety coated glass, 2 L
15 ml/s
2 ml - 500 ml
T-flasks, falcon tubes, petri dishes, multiwell culture plates, etc.
Yes (optional)

Glass, 0.5 L
2.3 ml/s
Up to 10 ml
Microcentrifuge tubes, microplates, small falcon tubes, small petri dishes, etc.
No

Part No.	Item Description	List Price	1 Unit 10% Discount	2-3 Units 20% Discount	4+ Units 30% Discount
158300	VACUSAFE 2 L bottle (glass) with bottle handle and closed lid Bottle adapter Level detection and lid with quick tubing couplings (PVDF) Tubing set with filter VACUBOY hand operator (155510)	3,250 CAD	2,925 CAD	2,600 CAD	2,275 CAD
158310	VACUSAFE Polypropylene bottle 4 L with lid Lid for Polypropylene bottle with quick couplings and level detection Tubing set with filter Mains adapter: 100-240VAC, 50/60Hz VACUBOY hand operator (155510)	2,600 CAD	2,340 CAD	2,080 CAD	1,820 CAD
158320	VACUSAFE Polypropylene bottle 4 L with lid Lid for Polypropylene bottle with barbed tube fittings Tubing set with filter Mains adapter: 100-240VAC, 50/60Hz VACUBOY hand operator (155510)	1,950 CAD	1,755 CAD	1,560 CAD	1,365 CAD
159000	VACUSIP Aspiration System with battery	895 CAD	806 CAD	716 CAD	627 CAD
159010	VACUSIP Aspiration System without battery	755 CAD	680 CAD	604 CAD	529 CAD



To benefit from this offer, please contact INTEGRA Biosciences by email or phone and mention the promotional code **Vacuum2024**, or clearly write the promotion code on faxed purchase orders. Discounts cannot be combined with any other offers. This offer is valid from May 1, 2024 through July 31, 2024.

INTEGRA Biosciences, US
22 Friars Drive
Hudson, NH 03051

Phone: 1-603-578-5800
Fax: 1-603-577-5529
E-Mail: us@integra-biosciences.com

www.integra-biosciences.com



Offer valid in Canada only.