

Improve the safety of your laboratory and save 30%!

VACUSAFE and **VACUSIP** liquid waste vacuum systems are the perfect replacement to in-house and self-made vacuum systems.

Both systems combine a vacuum pump, collection bottle, overflow protection and filter in an all-in-one liquid aspiration system. These systems offer the most effective method for containing biohazardous liquids with optimal personal protection.



Which system is right for me?



VACUSAFE



VACUSIP

	VACUSAFE	VACUSIP
Bottle Size	Polypropylene, 4 L/safety coated glass, 2 L	Glass, 0.5 L
Aspiration Range	15 ml/s	2.3 ml/s
Volumes for Aspiration	2 ml - 500 ml	Up to 10 ml
Typical Labware	T-flasks, falcon tubes, petri dishes, multiwell culture plates, etc.	Microcentrifuge tubes, microplates, small falcon tubes, small petri dishes, etc.
Level Detection (stops aspiration before overfilling)	Yes (optional)	No

Part No.	Item Description	List Price	30% Discount
158300	VACUSAFE 2 L bottle (glass) with bottle handle and closed lid Bottle adapter Level detection and lid with quick tubing couplings (PVDF) Tubing set with filter VACUBOY hand operator (155510)	CHF 2,830	CHF 1,981
158310	VACUSAFE Polypropylene bottle 4 L with lid Lid for Polypropylene bottle with quick couplings and level detection Tubing set with filter Mains adapter: 100-240VAC, 50/60Hz VACUBOY hand operator (155510)	CHF 2,338	CHF 1,637
158320	VACUSAFE Polypropylene bottle 4 L with lid Lid for Polypropylene bottle with barbed tube fittings Tubing set with filter Mains adapter: 100-240VAC, 50/60Hz VACUBOY hand operator (155510)	CHF 1,758	CHF 1,231
159000	VACUSIP Aspiration System with battery	CHF 813	CHF 569
159010	VACUSIP Aspiration System without battery	CHF 699	CHF 489



To benefit from this offer, please contact INTEGRABiosciences by email or phone and mention the promotional code **Vacuum2024**, or clearly write the promotion code on faxed purchase orders. Discounts cannot be combined with any other offers. This offer is valid from May 1, 2024 through July 31, 2024.

INTEGRA Biosciences AG
7205 Zizers, Switzerland

Phone: +41 81 286 95 55
Fax: +41 81 286 95 07
E-Mail: info-ch@integra-biosciences.com

www.integra-biosciences.com



Offer valid in Switzerland only.