

# Improve the safety of your laboratory and save up to 30%!

**VACUSAFE** and **VACUSIP** liquid waste vacuum systems are the perfect replacement to in-house and self-made vacuum systems.

Both systems combine a vacuum pump, collection bottle, overflow protection and filter in an all-in-one liquid aspiration system. These systems offer the most effective method for containing biohazardous liquids with optimal personal protection.



## Which system is right for me?



**VACUSAFE**



**VACUSIP**

	VACUSAFE	VACUSIP
<b>Bottle Size</b>	Polypropylene, 4 L/safety coated glass, 2 L	Glass, 0.5 L
<b>Aspiration Range</b>	15 ml/s	2.3 ml/s
<b>Volumes for Aspiration</b>	2 ml - 500 ml	Up to 10 ml
<b>Typical Labware</b>	T-flasks, falcon tubes, petri dishes, multiwell culture plates, etc.	Microcentrifuge tubes, microplates, small falcon tubes, small petri dishes, etc.
<b>Level Detection (stops aspiration before overfilling)</b>	Yes (optional)	No

Part No.	Item Description	List Price	1 Unit 10% Discount	2-3 Units 20% Discount	4+ Units 30% Discount
158300	<b>VACUSAFE</b> 2 L bottle (glass) with bottle handle and closed lid Bottle adapter Level detection and lid with quick tubing couplings (PVDF) Tubing set with filter VACUBOY hand operator (155510)	14.300,- DKK	12.870,- DKK	11.440,- DKK	10.010,- DKK
158310	<b>VACUSAFE</b> Polypropylene bottle 4 L with lid Lid for Polypropylene bottle with quick couplings and level detection Tubing set with filter Mains adapter: 100-240VAC, 50/60Hz VACUBOY hand operator (155510)	13.200,- DKK	11.880,- DKK	10.560,- DKK	9.240,- DKK
158320	<b>VACUSAFE</b> Polypropylene bottle 4 L with lid Lid for Polypropylene bottle with barbed tube fittings Tubing set with filter Mains adapter: 100-240VAC, 50/60Hz VACUBOY hand operator (155510)	10.000,- DKK	9.000,- DKK	8.000,- DKK	7.000,- DKK
159000	<b>VACUSIP</b> Aspiration System with battery	4.700,- DKK	4.230,- DKK	3.760,- DKK	3.290,- DKK
159010	<b>VACUSIP</b> Aspiration System without battery	3.900,- DKK	3.510,- DKK	3.120,- DKK	2.730,- DKK



To benefit from this offer, please contact INTEGRA Biosciences by email or phone and mention the promotional code **Vacuum2024**, or clearly write the promotion code on faxed purchase orders. Discounts cannot be combined with any other offers. This offer is valid from May 1, 2024 through July 31, 2024.

INTEGRA Biosciences Nordic ApS  
Vallensbækvej 22A 3TV  
2605 Brøndby  
Denmark

T +45 31 73 53 73

info-nordic@integra-biosciences.com  
www.integra-biosciences.com



Offer valid in Denmark only.