

Improve the safety of your laboratory and save up to 30%!

VACUSAFE and **VACUSIP** liquid waste vacuum systems are the perfect replacement to in-house and self-made vacuum systems.

Both systems combine a vacuum pump, collection bottle, overflow protection and filter in an all-in-one liquid aspiration system. These systems offer the most effective method for containing biohazardous liquids with optimal personal protection.



SAVE UP TO 30%!

Which system is right for me?



VACUSAFE

Bottle Size
Aspiration Range
Volumes for Aspiration
Typical Labware
Level Detection (stops aspiration before overfilling)

Polypropylene, 4 L/safety coated glass, 2 L
15 ml/s
2 ml - 500 ml
T-flasks, falcon tubes, petri dishes, multiwell culture plates, etc.
Yes (optional)



VACUSIP

Glass, 0.5 L
2.3 ml/s
Up to 10 ml
Microcentrifuge tubes, microplates, small falcon tubes, small petri dishes, etc.
No

Part No.	Item Description	List Price	1 Unit 10% Discount	2-3 Units 20% Discount	4+ Units 30% Discount
158300	VACUSAFE 2 L bottle (glass) with bottle handle and closed lid Bottle adapter Level detection and lid with quick tubing couplings (PVDF) Tubing set with filter VACUBOY hand operator (155510)	1,887 €	1,698 €	1,510 €	1,321 €
158310	VACUSAFE Polypropylene bottle 4 L with lid Lid for Polypropylene bottle with quick couplings and level detection Tubing set with filter Mains adapter: 100-240VAC, 50/60Hz VACUBOY hand operator (155510)	1,697 €	1,527 €	1,358 €	1,188 €
158320	VACUSAFE Polypropylene bottle 4 L with lid Lid for Polypropylene bottle with barbed tube fittings Tubing set with filter Mains adapter: 100-240VAC, 50/60Hz VACUBOY hand operator (155510)	1,267 €	1,140 €	1,014 €	887 €
159000	VACUSIP Aspiration System with battery	595 €	536 €	476 €	417 €
159010	VACUSIP Aspiration System without battery	495 €	446 €	396 €	347 €



To benefit from this offer, please contact INTEGRA Biosciences by email or phone and mention the promotional code **Vacuum2024**, or clearly write the promotion code on faxed purchase orders. Discounts cannot be combined with any other offers. This offer is valid from May 1, 2024 through July 31, 2024.

INTEGRA Biosciences SAS
Zone Industrielle des Béthunes - La Mare II
8 avenue du Fief
CS 79570 St Ouen l'Aumône
95062 Cergy Pontoise Cedex 1

Téléphone : +33 1 34 30 76 76
E-Mail : info-fr@integra-biosciences.com
www.integra-biosciences.com

INTEGRA

Offer valid in France only.