

Improve the safety of your laboratory and save 30%!

VACUSAFE and **VACUSIP** liquid waste vacuum systems are the perfect replacement to in-house and self-made vacuum systems.

Both systems combine a vacuum pump, collection bottle, overflow protection and filter in an all-in-one liquid aspiration system. These systems offer the most effective method for containing biohazardous liquids with optimal personal protection.



Which system is right for me?



VACUSAFE



VACUSIP

	VACUSAFE	VACUSIP
Bottle Size	Polypropylene, 4 L/safety coated glass, 2 L	Glass, 0.5 L
Aspiration Range	15 ml/s	2.3 ml/s
Volumes for Aspiration	2 ml - 500 ml	Up to 10 ml
Typical Labware	T-flasks, falcon tubes, petri dishes, multiwell culture plates, etc.	Microcentrifuge tubes, microplates, small falcon tubes, small petri dishes, etc.
Level Detection (stops aspiration before overfilling)	Yes (optional)	No

Part No.	Item Description	List Price	30% Discount
158300	VACUSAFE 2 L bottle (glass) with bottle handle and closed lid Bottle adapter Level detection and lid with quick tubing couplings (PVDF) Tubing set with filter VACUBOY hand operator (155510)	£1,813	£1,269
158310	VACUSAFE Polypropylene bottle 4 L with lid Lid for Polypropylene bottle with quick couplings and level detection Tubing set with filter Mains adapter: 100-240VAC, 50/60Hz VACUBOY hand operator (155510)	£1,671	£1,170
158320	VACUSAFE Polypropylene bottle 4 L with lid Lid for Polypropylene bottle with barbed tube fittings Tubing set with filter Mains adapter: 100-240VAC, 50/60Hz VACUBOY hand operator (155510)	£1,224	£857
159000	VACUSIP Aspiration System with battery	£540	£378
159010	VACUSIP Aspiration System without battery	£448	£314



To benefit from this offer, please contact INTEGRA Biosciences by email or phone and mention the promotional code **Vacuum2024**, or clearly write the promotion code on faxed purchase orders. Discounts cannot be combined with any other offers. This offer is valid from May 1, 2024 through July 31, 2024.

INTEGRA Biosciences Ltd.
2 Rivermead Business Park
Thatcham, Berks RG19 4EP, UK

Phone: +44 1635 797000
Fax: +44 1635 797001
E-Mail: info-uk@integra-biosciences.com

www.integra-biosciences.com



Offer valid in United Kingdom only.