

# Safe Disposal of Biological Liquid Waste

VACUSAFE Liquid Aspiration System and VACUBOY Vacuum Hand Operator

Save  
30 %!



## SAFE

- Hydrophobic filter protects the vacuum source
- Easy closing precision connectors to avoid the escape of aerosols
- Level sensor detects when bottle is full and switches off pump to prevent overfilling
- Stable stand to avoid spillage



## USER FRIENDLY

- Simply switch on, set the desired vacuum and start aspirating!
- Low noise level – below 65dB(A)
- Quick disconnect tubing enables an easy exchange of tubing sets



## VERSATILE

- VACUBOY adapter allows for aspiration of liquids from virtually any lab vessel!
- Adjustable vacuum pressure

Part No.	Product Description	List Price (1 Unit)	Discount Price 30 % Off
158320	VACUSAFE <ul style="list-style-type: none"> <li>• Polypropylene bottle, 4 liter with lid</li> <li>• Lid for polypropylene bottle with barbed tube fittings</li> <li>• VACUBOY hand operator (155510)</li> </ul>	1,661 CAD	<b>1,165 CAD</b>
158310	VACUSAFE <ul style="list-style-type: none"> <li>• Polypropylene bottle, 4 liter with lid</li> <li>• Lid for polypropylene bottle with quick couplings and level detection</li> <li>• VACUBOY hand operator (155510)</li> </ul>	2,207 CAD	<b>1,545 CAD</b>

\*Additional accessories and VACUBOY hand operator adaptors available at retail rate

# Safe Disposal of Biological Liquid Waste

## VACUSIP Portable Liquid Aspiration System

Save  
30 %!



### SAFE

- Hydrophobic filter protects the vacuum source
- Easy decontamination via autoclave



### PORTABLE

- Compact, lightweight design
- Fits in nearly all biosafety cabinets
- Rechargeable battery



### CONVENIENT

- Minimal tubing/cables
- Vacuum pump runs silently
- Integrated holder for handheld aspiration device

Part No.	Product Description	List Price (1 Unit)	Discount Price 30 % Off
159000	VACUSIP Aspiration System with battery	824 CAD	575 CAD
159010	VACUSIP Aspiration System without battery	690 CAD	485 CAD

\*Additional accessories and VACUBOY hand operator adaptors available at retail rate



Check out this article on how the VACUSIP is being used with COVID-19 testing!

To benefit from this offer, please contact INTEGRA Biosciences by email or phone and mention the promotion code **VACU2020**, or clearly write the promotion code on all faxed purchase orders. This offer is valid until **December 31, 2020**.

Offer valid in Canada only.

INTEGRA Biosciences, US  
22 Friars Drive  
Hudson, NH 03051

Phone: 1-603-578-5800  
Fax: 1-603-577-5529  
E-Mail: [us@integra-biosciences.com](mailto:us@integra-biosciences.com)

[www.integra-biosciences.com](http://www.integra-biosciences.com)

INTEGRA