



PIPETBOY acu 2 – Pipette Controllers

The Most Popular Pipette Controller

Part No.	Description	List Price	Promo Price
155000	PIPETBOY acu 2 classic, incl. wall mount, mains adapter, sterile filter, battery	\$312	\$219
155015	PIPETBOY acu 2 green, as item no. 155 000	\$312	\$219
155016	PIPETBOY acu 2 red, as item no. 155 000	\$312	\$219
155017	PIPETBOY acu 2 blue, as item no. 155 000	\$312	\$219
155018	PIPETBOY acu 2 transparent, as item no. 155 000	\$312	\$219
155019	PIPETBOY acu 2 purple, as item no. 155 000	\$312	\$219
155022	PIPETBOY acu 2 Agave	\$312	\$219
155023	PIPETBOY acu 2 Ocean Dream	\$312	\$219
155024	PIPETBOY acu 2 Pink Sunrise	\$312	\$219



FIREBOY – Safety Bunsen Burners

A Safer Alternative to the Traditional Bunsen Burner


Part No.	Description	List Price	Promo Price
144000	FIREBOY plus – safety Bunsen burner with display, movement sensor, flame monitor and overheating protection. three operation modes (manual sensor, foot switch or continuous). Includes: <ul style="list-style-type: none"> • 2 nozzles for propane/butane (P60) and natural gas (N80) • Mains adapter • Adapter for standard gas hose with inner diameter 10 mm 	\$811	\$568
144010	FIREBOY eco – safety Bunsen burner with foot switch, including: <ul style="list-style-type: none"> • Mains adapter • 2 nozzles for propane/butane (P60) and natural gas (N80) • Adapter for standard gas hose with inner diameter 10 mm 	\$705	\$494



FLAMEBOY

Surface Sterilization at the Push of a Button

Part No.	Description	List Price	Promo Price
145000	FLAMEBOY – automatic flame pistol for propane, butane and town gas. Includes: <ul style="list-style-type: none"> • 3 nozzles for the various types of gas • Adapter for standard gas hose with inner diameter 10 mm • Threaded adapter for 1/4" gas hose, left hand thread 	\$525	\$368
145003	FLAMEBOY – automatic flame pistol for propane, butane and town gas with adapter for gas cartridge. Includes: <ul style="list-style-type: none"> • 3 nozzles for the various types of gas • Adapter for gas cartridge Campinggaz CV 360 	\$650	\$455



VACUBOY

Convert any Vacuum Source into a User-Friendly Aspiration System

30%
OFF

Part No.	Description	List Price	Promo Price
155500	VACUBOY hand operator – set consists of: <ul style="list-style-type: none"> VACUBOY hand operator Stand for hand operator Rubber adapter for Pasteur pipettes 1 channel stainless steel aspiration adapter 40 mm 8 channel stainless steel aspiration adapter 40 mm 1 channel plastic aspiration adapter for disposable tips 1 channel plastic aspiration adapter for disposable tips with ejector 8-channel plastic aspiration adapter for disposable tips with ejector 	\$375	\$263
155510	VACUBOY hand operator with 1 channel stainless aspiration adapter 40 mm and 1 channel plastic aspiration adapter for disposable tips with ejector, incl. rubber adapter for Pasteur pipettes.	\$198	\$139



30%
OFF

DOSE IT

An Easy-to-Use Peristaltic Pump

Part No.	Description	List Price	Promo Price
171000	DOSE IT P910 programmable peristaltic pump, including: <ul style="list-style-type: none"> Silicone tubing 4 mm inner diameter (ID) 2 Aspiration/dispensing tubes, length 10 cm, stainless steel, one end dented 1 Aspiration/dispensing tubes, length 35 cm, stainless steel, one end dented Tube collar 1 Pipette-tubing connector Mains adapter 	\$2500	\$1750



30%
OFF

EVOLVE – Manual Pipettes

Set volumes in the blink of an eye instead of a twist of the wrist!



Part No.	Description	List Price	Promo Price
3011	Single channel, 0.2-2 µl	\$355	\$249
3012	Single channel, 1-10 µl		
3013	Single channel, 2-20 µl		
3015	Single channel, 10-100 µl		
3016	Single channel, 20-200 µl	\$335	\$235
3018	Single channel, 100-1000 µl		
3019	Single channel, 500-5000 µl		
3022	8 channel, 1-10 µl	\$760	\$532
3024	8 channel, 5-50 µl		
3025	8 channel, 10-100 µl	\$740	\$518
3026	8 channel, 20-200 µl		

Part No.	Description	List Price	Promo Price
3028	8 channel, 100-1250 µl	\$835	\$585
3032	12 channel, 1-10 µl	\$865	\$606
3034	12 channel, 5-50 µl		
3035	12 channel, 10-100 µl	\$845	\$592
3036	12 channel, 20-200 µl		
3038	12 channel, 100-1250 µl	\$950	\$665
3042	16 channel, 1-10 µl		
3044	16 channel, 5-50 µl	\$995	\$697
3045	16 channel, 10-100 µl		

To benefit from this offer, please contact INTEGRA Biosciences by email or phone and mention the promotion code **USEndofYear**, or clearly write the promotion code on all faxed purchase orders. This offer is valid until **December 31st 2020**.

Offer valid in the United States only.

INTEGRA Biosciences, US
2 Wentworth Drive
Hudson, NH 03051

Phone: 1-603-578-5800
Fax: 1-603-577-5529
E-Mail: us@integra-biosciences.com

www.integra-biosciences.com

INTEGRA