


INTEGRA

“The Bead has landed”



Clean Blood & Tissue DNA Kit Magnetic bead-based genomic DNA extraction kit

Houston, we have a solution.
DNA extraction without centrifugation, so it's ready to automate.

Get your free kit now!

 bit.ly/cbtkit-sample



A better way to extract genomic DNA

Is your centrifuge busy? The Clean Blood & Tissue DNA Kit uses magnetic beads instead of spin columns, so you don't need a centrifuge.


Extract genomic DNA with a workflow that's simple enough for manual processing and ready for automation.

The workflow

Lyse > Bind > Wash > Elute

Cells are lysed to release genomic DNA, which binds to magnetic beads. Magnets hold the beads in place during washing, then the purified DNA is released during elution.

Automation when you need it



Ready to automate? Use the ASSIST PLUS for hands-free extraction with improved reproducibility across runs and users.

Ordering information

Product	Preps	Part no.
Clean Blood & Tissue DNA Kit	16	7100
	96	7101
	384	7102

Your key benefits



Consistently high yields

from blood, buffy coat, tissue and cultured cells



High purity DNA

ready for (q)PCR and NGS



Automation ready

to save hands-on time and improve reproducibility

Explore the kit!

 bit.ly/cbt-kit



INTEGRA Biosciences AG
Tardisstrasse 201
CH-7205 Zizers
Switzerland

Phone: +41 81 286 95 55
E-Mail: info-ch@integra-biosciences.com

www.integra-biosciences.com



INTEGRA