

### EVOLVE – Manual Pipettes

Adjust volumes more than ten times faster



Single Channel	
0.2	- 2 $\mu$ l
1	- 10 $\mu$ l
2	- 20 $\mu$ l
10	- 100 $\mu$ l
20	- 200 $\mu$ l
100	- 1,000 $\mu$ l
500	- 5,000 $\mu$ l

#### 8, 12 and 16 Channel

1	- 10 $\mu$ l (8, 12, 16)
5	- 50 $\mu$ l (8, 12, 16)
10	- 100 $\mu$ l (8, 12, 16)
20	- 200 $\mu$ l (8, 12)
100	- 1,250 $\mu$ l (8, 12)



### VIAFLO – Electronic Pipettes

GRIPTIPS® pipette tips will never loosen, leak or fall off



Single Channel	
0.5	- 12.5 $\mu$ l
2	- 50 $\mu$ l
5	- 125 $\mu$ l
10	- 300 $\mu$ l
50	- 1,250 $\mu$ l
100	- 5,000 $\mu$ l


#### 8, 12 and 16 Channel

0.5	- 12.5 $\mu$ l (8, 12, 16)
2	- 50 $\mu$ l (8, 12, 16)
5	- 125 $\mu$ l (8, 12, 16)
10	- 300 $\mu$ l (8, 12)
50	- 1,250 $\mu$ l (8, 12)



### VOYAGER – Electronic Tip Spacing Pipettes

Transferring samples between different labware



#### 4, 6, 8 and 12 Channel

0.5	- 12.5 $\mu$ l (8, 12)
2	- 50 $\mu$ l (8, 12)
5	- 125 $\mu$ l (8, 12)
10	- 300 $\mu$ l (4, 6, 8)
50	- 1,250 $\mu$ l (4, 6, 8)



### MINI 96, VIAFLO 96, and VIAFLO 384 – Electronic Pipettes

Significant increase in pipetting productivity, speed and reproducibility



#### MINI 96 96 channels

0.5	- 12.5 $\mu$ l
5	- 125 $\mu$ l
10	- 300 $\mu$ l
50	- 1,250 $\mu$ l

#### VIAFLO 96 and VIAFLO 384 24, 96, and 384 channels

0.5	- 12.5 $\mu$ l (96, 384)
2	- 50 $\mu$ l (96, 384)
5	- 125 $\mu$ l (96, 384)
10	- 300 $\mu$ l (24, 96)
50	- 1,250 $\mu$ l (24, 96)



### ASSIST and ASSIST PLUS – Pipetting Robots

Automate INTEGRA's VIAFLO and VOYAGER pipettes and the D-ONE single channel pipetting modules\*

#### Applications:

- Serial dilutions
- Plate fillings
- Reagent additions
- Sample reformatting\*
- Hit picking\*
- Normalization\*
- Master mix\*



ASSIST



D-ONE

Single channel pipetting module for ASSIST PLUS



ASSIST PLUS

\* Only available for ASSIST PLUS

## MAG, HEATMAG, COLDPLATE and BIOSHAKE

Magnetic separation, heating, cooling and shaking modules



**MAG**

Automated magnetic bead collection and capture



**HEATMAG**


Combined automatic magnetic bead collection and heating function

Compatible with INTEGRA's benchtop pipetting platforms



**COLDPLATE**

Precise and uniform temperature control from 0 to 99 °C



**BIOSHAKE**

Precise shaking up to 3000 rpm and accurate heating



## WELLJET – Reagent Dispenser And Dispenser Stacker

Low cost, fast, and precise bulk reagent dispensing



**Dispenser**

For rapid filling of 6 to 1536 well plates  
0.5 µl – 10 ml



**Dispenser Stacker**

For hands-free filling of non lidded and lidded well plates

## PIPETBOY acu 2, PIPETBOY pro – Pipet Controller

Precise and efficient serological pipetting

**PIPETBOY acu 2**

Fits serological pipettes from 1 ml - 100 ml



**PIPETBOY pro**

With LED light



## VACUSAFE, VACUBOY, VACUSIP, DOSE IT

Liquid aspiration, collection, disposal and pump systems



**VACUSAFE and VACUBOY**

Safe aspiration system with 2 L and 4 L bottles

**VACUSIP**

Portable aspiration system with 500 ml bottle



**DOSE IT**

Easy to operate peristaltic pump



## MIRO CANVAS – NGS Prep System

Enables full automation of NGS prep protocols



**MIRO CANVAS**

Compact microfluidics platform that enables true walk-away automation of complex NGS library prep and hybridization protocols

**MIRO CARTRIDGE**

All protocol steps, including mixing, thermal cycling and magnetic bead clean-ups are performed in the MIRO cartridge



## Eco-Friendly Reservoirs, Smart Tips

Everything you need for productive pipetting from 1 to 384 channels

Wide range of reservoirs with reusable bases and disposable inserts



GRIP TIPS® never loosen, leak or fall off