

Product name: GRIPTIPS 300 µl long
AUTO rack of 96 tip sterile
Lot Number: all lots

Article number: 6484

Our commitment to customers:

INTEGRA Biosciences polymer products are using precise and consistent manufacturing techniques that ensure quality and excellent performance. We strive to exceed our customers' expectations and encourage continuous feedback. The reservoirs are manufactured in accordance to ISO 9001.

Our product lots are traceable down to the corresponding raw materials used and date produced.

Our manufacturing process incorporates consistent and continuous inspection throughout all phases of the manufacturing process. The reservoirs produced under these strict conditions are contaminant free and consistent from lot-to-lot.

High Sensitivity Endotoxin Testing:

For products used in endotoxin testing we have performed additional testing for the presence of endotoxins.

Our Swiss reference lab, accredited to ISO 17025, performed additional lot-specific Limulus amoebocyte lysate (LAL) test using the turbidimetric test procedure conforming to EP 2.6.14 / USP <85>. The cut-off level for endotoxin acceptance is set at at: 0.005 EU/ml.



Urs Hartmann
CEO
INTEGRA Biosciences AG



Dr. Daniel Bächli
Head of Corporate Quality
INTEGRA Biosciences AG



INTEGRA Biosciences AG
7205 Zizers, Switzerland
T +41 81 286 95 30
info@integra-biosciences.com

INTEGRA Biosciences Corp.
Hudson NH 03051, USA
T +1 603 578 5800
info-us@integra-biosciences.com

INTEGRA Biosciences Deutschland GmbH
35444 Biebertal, Deutschland
T +49 6409 81 999 15
info-de@integra-biosciences.com

INTEGRA Biosciences SAS
95062 Cergy-Pontoise Cedex 1, France
T +33 1 34 30 76 76
info-fr@integra-biosciences.com

INTEGRA Biosciences Ltd.
Thatcham Berkshire RG194EP, UK
T +44 1635 797 00
info-uk@integra-biosciences.com

INTEGRA Biosciences (Shanghai) Co., Ltd.
上海自由贸易试验区 201315, 中国
T +86 21 5844 7203
info-cn@integra-biosciences.com

INTEGRA Biosciences KK
ちよだ一く、ときよ 101-0031、じゃぼん
T +81 3 5962 4936
info-jp@integra-biosciences.com

INTEGRA Biosciences Nordic ApS
2605 Brøndby, Denmark
T +45 3173 5373
orders-nordic@integra-biosciences.com