














## Overview: Supported inkjet printers (1)

| Supported printers:  | MI 9330  | MI 9029  | MI 9410   | MI 9450   |
|--|--|--|---|---|
| <ul style="list-style-type: none"> <li>MI 9029</li> <li>MI 9330</li> <li>MI 9410</li> <li>MI 9450</li> <li>MI 9018 (discont.)</li> <li>MI 9020 (discont.)</li> <li>MI 9028 (discont.)</li> <li>MI 9030 (discont.)</li> <li>MI 9232 (discont.)</li> <li>MI 9341 (discont.)</li> </ul> | <br><a href="http://www.markem-ijaje.com/">www.markem-ijaje.com/</a>                          | <br><a href="http://www.markem-ijaje.com/">www.markem-ijaje.com/</a>                         | <br><a href="http://www.markem-ijaje.com/">www.markem-ijaje.com/</a>                           | <br><a href="http://www.markem-ijaje.com/">www.markem-ijaje.com/</a>                           |
| <b>Performance</b>   | <ul style="list-style-type: none"> <li>4 Lines</li> <li>Alphanumeric</li> <li>wide choice of characters</li> <li>Graphics/Logo printing</li> <li>Side wall printing</li> </ul> | <ul style="list-style-type: none"> <li>4 Lines</li> <li>Alphanumeric</li> <li>wide choice of characters</li> <li>Graphics/Logo printing</li> <li>Side wall printing</li> </ul> | <ul style="list-style-type: none"> <li>5 Lines</li> <li>Alphanumeric, bar codes, 2D codes</li> <li>Graphics/Logo printing</li> <li>Side wall printing</li> </ul>                  | <ul style="list-style-type: none"> <li>5 Lines</li> <li>Alphanumeric, bar codes, 2D codes</li> <li>Graphics/Logo printing</li> <li>Side wall printing</li> </ul>                  |
| <b>Other features</b>  | <ul style="list-style-type: none"> <li>Bench top</li> <li>Touch screen</li> <li>USB port</li> <li>RS-232 / Ethernet</li> </ul>   | <ul style="list-style-type: none"> <li>Bench top</li> <li>Key pad</li> <li>USB port</li> <li>RS-232</li> </ul>   | <ul style="list-style-type: none"> <li>Bench top</li> <li>Touch screen</li> <li>USB and SD port</li> <li>Ethernet</li> </ul>  | <ul style="list-style-type: none"> <li>Bench top</li> <li>Touch screen</li> <li>USB and SD port</li> <li>Ethernet</li> </ul>  |
| <b>Dimensions (HxWxD)</b>  | 595 x 303 x 280 mm   | 595 x 307 x ??? mm   | 555 x 360 x 357 mm  | 555 x 360 x 357 mm  |
| <b>Weight (kg)</b>   | 18.5   | 18.5   | 25  | 25  |
| <b>List Price (EUR)</b>  | <10'000  | <10'000  | >10'000   | >10'000   |
| <b>Maintenance</b>   | <ul style="list-style-type: none"> <li>ink exchange in &lt; 6 min</li> <li>ink runs &gt; 8000 hours</li> </ul>   | <ul style="list-style-type: none"> <li>ink exchange in &lt; 6 min</li> <li>ink runs &gt; 8000 hours</li> <li>automatic cleaning system</li> </ul>                              | <ul style="list-style-type: none"> <li>Maintenance: every 18 month</li> <li>Head cleaning recommended every second week</li> <li>Head can be left attached on MEDIAJET</li> </ul> | <ul style="list-style-type: none"> <li>Maintenance: every 18 month</li> <li>Head cleaning recommended every second week</li> <li>Head can be left attached on MEDIAJET</li> </ul> |
| <b>INTEGRA comments</b>  | <b>Recommended MI printer</b> , lower consumable and maintenance costs. Industrial interface must be ordered.  | Industrial interface must be ordered.  | Industrial interface must be ordered.   | Industrial interface must be ordered. Higher line speed and storage capacity.   |

## Overview: Supported inkjet printers (2)

| Supported printers:  | Linx 89xx<br>(89xx: -00, -10, -20, -40)   | Domino Ax150i  | Domino Ax350i  |
|--|---|--|--|
| <ul style="list-style-type: none"> <li>• Linx 4900, Linx 5900 (printhead MK7)</li> <li>• Linx CJ400, Linx 89xx (printhead MK11)</li> <li>• Domino A320i, Domino A420i, Domino Ax150i, Domino Ax350i</li> </ul> | <br><a href="http://www.linxglobal.com">www.linxglobal.com</a>                               | <br><a href="http://www.domino-printing.com">www.domino-printing.com</a>                       | <br><a href="http://www.domino-printing.com">www.domino-printing.com</a>                      |
| Performance  | <ul style="list-style-type: none"> <li>▪ 6 Lines</li> <li>▪ Alphanumeric, barcodes, 2D codes</li> <li>▪ Graphics/Logo 25 px</li> <li>▪ Side wall / bottom printing</li> </ul> | <ul style="list-style-type: none"> <li>▪ 3 Lines</li> <li>▪ Alphanumeric, bar codes, 2D codes</li> <li>▪ Graphics/Logo printing</li> <li>▪ Side wall printing</li> </ul>         | <ul style="list-style-type: none"> <li>▪ 5 Lines</li> <li>▪ Alphanumeric, bar codes, 2D codes</li> <li>▪ Graphics/Logo printing</li> <li>▪ Side wall printing</li> </ul>         |
| Other features   | <ul style="list-style-type: none"> <li>▪ Benchtop</li> <li>▪ Touch screen</li> <li>▪ USB port</li> <li>▪ Ethernet</li> </ul>  | <ul style="list-style-type: none"> <li>▪ Bench top</li> <li>▪ Touch screen</li> <li>▪ Optional: USB port, Ethernet, RS232, GPIO</li> </ul>                                       | <ul style="list-style-type: none"> <li>▪ Bench top</li> <li>▪ Touch screen</li> <li>▪ Optional: USB port, Ethernet, RS232, GPIO</li> </ul>                                       |
| Dimensions (HxWxD)   | 503 x 440 x 355 mm  | 373 x 400 x 393 mm   | 411 x 430 x 331 mm   |
| Weight (kg)  | 24  | 16   | 22   |
| List Price (EUR)   | >10'000   | 8'300  | >10'000  |
| Maintenance  | <ul style="list-style-type: none"> <li>▪ Maintenance: every 24 months</li> <li>▪ Head cleaning every 3 months</li> </ul>  | <ul style="list-style-type: none"> <li>▪ after 3500 working hours no external service necessary, 10 min. maintenance</li> <li>▪ Head can be left attached on MEDIAJET</li> </ul> | <ul style="list-style-type: none"> <li>▪ after 3500 working hours no external service necessary, 10 min. maintenance</li> <li>▪ Head can be left attached on MEDIAJET</li> </ul> |
| INTEGRA comments   | This reliable Linx printer replaces all previous models.  | <b>Recommended Domino printer</b><br>Recommended, if bar codes or 2D codes needed.   | Recommended, if bar codes or 2D codes needed.  |

## Ordering information MEDIAJET inkjet printer connections kits

| Supported printers:  | Markem-Imaje (MI)<br>9018-9030, 9330  | Markem-Imaje (MI)<br>9232/94xx  | Linx<br>side wall printing<br>(with print head MK7<br>and MK11)                     | Linx<br>bottom printing<br>(with print head MK7<br>and MK11)                        | Domino  |
|--|---|---|---|---|---|
| <ul style="list-style-type: none"> <li>MI 9018, MI 9020, MI 9028, MI 9029, MI 9030, MI9330</li> <li>MI 9232, MI 9410, MI 9450</li> <li>Linx 4900, Linx 5900, Linx 89xx (print head MK7)</li> <li>Linx CJ400 (print head MK11)</li> <li>Domino A320i, Domino A420i, Domino Ax150i, Domino Ax350i</li> </ul> |    |  |  |  |  |
| <b>Part number</b>   |   |   |   |   |   |
| <b>Printer connection kit:</b><br>with printer head fixing device and interface cable<br>   | 103080  | 103091  | 113845  | 113840  | 113060  |
| <b>Operating instructions</b>  | Operating instructions for connection kits are included and also available in the download center on <a href="http://www.integra-biosciences.com">www.integra-biosciences.com</a> |   |   |   |   |
|  | 103983  | 103980  | 103956  | 103981  | 103982  |
| <b>Cart</b><br>(W x D x H (cm):<br>70 x 50 x 65-90)<br>  | -   | -   | -   | 113841  | -   |
| <b>INTEGRA comments</b>  | Compatible with:<br>Ø 90 mm, Ø 60 mm or<br>Ø 35 mm Petri dishes.  | Compatible with:<br>Ø 90 mm or Ø 60 mm<br>Petri dishes.                             | Compatible with:<br>Ø 90 mm or Ø 60 mm<br>Petri dishes.                             | Compatible with:<br>Ø 90 mm, Ø 60 mm or<br>Ø 35 mm Petri dishes.                    | Compatible with:<br>Ø 90 mm or Ø 60 mm<br>Petri dishes.                             |
| <b>Compatibility with MEDIAJET (manufactured)</b>  | After Dec 2005 / Feb 2007; (SN above:<br>MJ : 10305000 /<br>MJ vario : 11307000)  | After Feb 2013;<br>(SN above:<br>MJ : 10313568 /<br>MJ vario : 11313247)            | After Feb 2013;<br>(SN above:<br>MJ : 10313568 /<br>MJ vario : 11313247)            | After Mai 2010;<br>(SN above:<br>MJ : 10310371 /<br>MJ vario : 11310128)            | After Feb 2013;<br>(SN above:<br>MJ : 10313568 /<br>MJ vario : 11313247)            |