

Reagent reservoirs: Divided reagent reservoirs

Laboratory equipment | Labware | Reagent Reservoirs



25 ml divided reservoir with lowest dead volume

Product description:

25 ml divided reservoirs are comprised of 5 and 10 ml compartments, while still nesting comfortably in the standard 25 ml reusable, sturdy reservoir base. Each reservoir features the SureFlo™ anti-sealing array and hydrophilic surface treatment providing especially low dead volumes. This is particularly helpful when pipetting precious samples or expensive reagents with a multichannel pipette. Inserts are available in polystyrene or polypropylene for improved chemical compatibility.

Features and benefits:

- Lowest dead volume due to SureFlo™ anti-sealing pattern, inert hydrophilic surface treatment and low volume compartments of 5 and 10 ml.
- Less plastic waste, divided inserts nest in the standard 25 ml reusable base.
- Save on reagents with pour-back spouts moulded into the corners to return excess reagent to its original container.
- Base accommodates two reservoirs, one acting as a lid for short term storage and to prevent contamination.