

Pipette Tips: ECO Rack GRIPTIPS for Handheld Pipettes

Laboratory equipment | Pipetting devices | Pipette tip



ECO Rack GRIPTIPS for the perfect connection between pipette and tip

Product description:

The GRIPTIPS together with INTEGRA's handheld pipettes form an integral system. Pipette tip and pipette were designed to fit perfectly together. As a result, the loading and attachment forces of the tips are reduced to minimize wrist related stresses. In addition, the system provides a perfect tip sealing and reliably prevents tips from accidentally loosening, leaking or falling off during pipetting.

Available in five volumes (12.5 μ l, 50/125 μ l, 300 μ l, 1250 μ l and 5000 μ l) GRIPTIPS come in bulk or stackable boxes (non-sterile, sterile and filtered, low retention) featuring color-coded inserts that correspond to the volume ranges of the INTEGRA pipettes. The 12.5 μ l and 125 μ l GRIPTIPS come in 384-fold configuration, allowing the same tips and boxes to be used with 16-channel VIAFLO pipettes.

The ECO Rack is an environmentally friendly packaging of GRIPTIPS with 60% less plastic than alternative racks. The reusable PopTop base is designed to securely hold ECO Racks and GREEN CHOICE GRIPTIPS for an optimal pipette tip loading experience.

Features and benefits:

- GRIPTIPS snap onto the tip fitting of INTEGRA's pipettes, which prevents tips from falling off.
- Reduced tip attachment and ejection forces for comfortable, stress-free pipetting.
- Precise and consistent tip seal ensures that all tips on a multichannel pipette are precisely at an identical height for reliable touch offs and consistent pipetting results.
- Reduce plastic waste without sacrificing quality.