

Pipette controller: PIPETBOY pro

Laboratory equipment | Pipetting devices | battery powered pipette controller



PIPETBOY pro – Enlightened pipetting

Product description:

Accuracy, speed and ergonomic design are the essential characteristics realized in the PIPETBOY pro in order to productively work with all serological pipettes of 1-100 ml.

Designed for error-free and comfortable pipetting, the white LED light integrated in the PIPETBOY pro provides an optimal illumination of the area around the pipette, ensuring accurate pipetting and preventing eyestrain.

A highly accurate flow control is achieved simply by applying variable finger pressure on the pipettor trigger buttons, changing between drop-by-drop dispensing or fast motorized pipetting at speeds of up to 10 ml/s. Weighing just 190 grams and ergonomically designed to comfortably fit in your hand, the PIPETBOY pro makes even prolonged pipetting effortless.

Features and benefits:

- Integrated LED illumination of pipette for error- and eyestrain-free working.
- Fine control of pipetting via finger pressure enables accurate drop-by-drop dispensing.
- Quiet and strong pump for high pipetting speed of up to 10 ml/s.
- Lightweight and ergonomic design for maximal comfort.
- Removable and autoclavable pipette holder module for highest cleanness and safety.
- Five colours to choose from according to personal preferences and lab organization.