

Aspiration systems: VACUBOY

Laboratory equipment | Aspirators



VACUBOY – For aspiration of liquids from all sorts of containers

Product description:

The VACUBOY is a convenient handheld tool that converts a vacuum source into a user friendly system for aspirating liquids from virtually any laboratory container – at the touch of a button and in a controlled way. The hand piece, the core element of the VACUBOY, is simply connected to any available vacuum source using silicone tubing and fitted with one of the different aspiration adapters. The perfect match is achieved when the VACUBOY is connected to a VACUSAFE aspiration system.

Features and benefits:

- Adapters for aspiration with Pasteur/serological pipettes, disposable pipette tips or stainless steel needles.
- Multi-channel adapters for fast emptying of multi-well plates.
- Completely autoclavable for maximal safety in the lab.
- Vacuum regulation by a touch-sensitive button for controlled and smooth aspiration of liquids.