

Electronic Pipettes: VIAFLO Electronic Pipettes

Laboratory equipment | Pipetting devices | Multichannel electronic pipette



VIAFLO – Lightweight Multichannel Electronic Pipettes

Product description:

The single, 8-, 12- and 16-channel models of VIAFLO electronic pipettes together with INTEGRA's wide range of GRIPTIPS form a perfect pipetting system. GRIPTIPS snap into place with minimal tip loading effort, providing a secure connection. GRIPTIPS never fall off and are always perfectly aligned, resulting in superior accuracy and precision results. The significantly reduced, well balanced weight of the VIAFLO pipettes features unsurpassed multichannel pipetting ergonomics, even during prolonged pipetting sessions. The color screen together with the easy to navigate touch wheel user interface promote ease of set up and operation. All VIAFLO pipettes offer a choice of predefined pipetting modes covering typical applications as well as custom multi-step programs for specific tasks. VIAFLO pipettes are available with volume ranges of 0.5 – 12.5 μl , 2 – 50 μl , 5 – 125 μl , 10 – 300 μl and 50 – 1250 μl and 500 – 5000 μl for the single channel VIAFLO model.

Features and Benefits:

- VIAFLO pipettes are up lightweight and help to reduce risk of RSI injuries.
- Robust and easy to repair design promotes trouble free pipetting.
- Additional menu languages in German and French, Spanish, Chinese and Japanese.
- Multichannel VIAFLO pipettes can be attached to ASSIST and ASSIST PLUS to perform pipetting operations automatically.