

Microplate dispensers: WELLJET

Laboratory equipment | Laboratory automation equipment | automatic microplate dispenser



WELLJET - Bulk reagent dispenser and dispenser stacker

Product description:

The WELLJET dispenser and dispenser stacker are designed to offer superior user experience and flexibility for applications requiring low cost and precise bulk reagent dispensing. Both models use a revolutionary type of dispensing cassettes (EasySnap™) which allow for low running costs and very easy handling. With their highly compact design, both the dispenser and the dispenser stacker exhibit the smallest footprints for instruments of their kind, saving valuable bench space. A large, high resolution touchscreen user interface guarantees intuitive and effortless programming.

Features and benefits:

- Compact footprint – both fit into a laminar flow cabinet
- Fast filling of microplates
- Low running costs with the EasySnap™ cassette
- More walkaway time with robust and reliable, hands-free filling of lidded and unlidded plates with the dispenser stacker