

| | WELLJET dispenser | WELLJET dispenser Stacker |
|---------------------------------|---|---|
| Temperature range | 5 - 35°C | |
| Humidity range | Max. rel. humidity 80% for temperatures up to 31°C, decreasing linearly to 50% rel. humidity at 40°C. | |
| Dimensions (W x D x H) | 20 x 46 x 29 cm (7.9 x 11.4 x 18.1 in) | 46 x 46 x 63/102* cm (18.1 x 18.1 x 24.8/40.2* in) * tower for 25/50 plates |
| Weight | 8.8 kg (19.4 lbs) | 22.1 kg (48.7 lbs) with 2 towers for 25 plates |
| Electricity supply | Mains adapter input: 100–240 VAC, 47–63 Hz Device input: 22.8–25.2 VDC, 100 W | |
| Noise emission | <60 dBA | |
| User interface | 17.8 cm / 7" high resolution touch screen | |
| Integration | USB interface API commands available | - |
| Detection method | - | optical sensor |
| Compatible plate formats | 6, 12, 24, 48, 96, 384, 1536, shallow & deep well plates | |
| Compatible plate heights | 5-64 mm | 9-25 mm (stacker mode) |
| | | 5-64 mm (manual mode) |