

Kanal	Art.-Nr.	Volumen-Bereich (µl)	Volumen-schritt (µl)	Test-Volumen (µl)	Genauig-keit (±%)	Präzision (≤%)
1	4531	0.5 - 20	0.01	2	10.00	6.00
				10	4.00	1.60
				20	2.00	0.80
	10 - 300	0.5	30	4.00	1.20	
			150	2.00	0.60	
			300	1.60	0.35	
	4532	5 - 125	0.1	12.5	3.75	1.50
				62.5	2.50	0.70
125				1.60	0.35	
50 - 1250	1	1	125	6.00	1.10	
			625	2.40	0.50	
			1250	1.20	0.30	

Präzision = Streukoeffizient