

Specifications of D-ONE pipetting modules



Channel	Part No.	Volume Range (µl)	Volume Increments (µl)	Test Volume (µl)	Accuracy (±%)	Precision (≤%)
1	4531	0.5 - 20	0.01	2	10.00	6.00
				10	4.00	1.60
				20	2.00	0.80
	4532	10 - 300	0.5	30	4.00	1.20
				150	2.00	0.60
				300	1.60	0.35
				5 - 125	0.1	12.5
62.5	2.50	0.70				
125	1.60	0.35				
50 - 1250	1	125	6.00	1.10		
		625	2.40	0.50		
		1250	1.20	0.30		

Precision = coefficient of variation