

Canaux	Réf. Art.	Gamme de volume (µl)	Incément de volume (µl)	Volume test (µl)	Exactitude (±%)	Précision (≤%)
1	4531	0.5 - 20	0.01	2	10.00	6.00
				10	4.00	1.60
				20	2.00	0.80
		10 - 300		30	4.00	1.20
				150	2.00	0.60
				300	1.60	0.35
	4532	5 - 125	0.1	12.5	3.75	1.50
				62.5	2.50	0.70
				125	1.60	0.35
		50 - 1250		1	125	6.00
625	2.40		0.50			
1250	1.20	0.30				

Précision = Coefficient de Variation