


Get your pipettes calibrated by INTEGRA – on site, directly in your lab!

INTEGRA strongly recommends an annual calibration to ensure peak performance of your pipettes. The best way for you to accomplish this is to have INTEGRA's certified calibration partner come to your lab and perform your pipette calibrations on site.

Guarantee that the performance of your pipettes meets your highest expectations!

- INTEGRA's on-site calibration process is much faster and easier than sending pipettes out to a calibration lab, and requires less paperwork and no packing
- Pipettes are calibrated in your own lab climate conditions
- Fastest possible turnaround time - your pipettes can be used immediately.



INTEGRA's on-site pipette calibration service

Our partner performs calibration of all of your manual and electronic pipettes, whether they are INTEGRA pipettes or from other major manufacturers. All you need to provide is the pipettes, and we bring the rest. This service guarantees that your pipettes will perform at their best, while being the easiest solution for you.


INTEGRA recommends an annual service frequency.

On-site calibration

Visual inspection for external damage and cleaning	✓
Change of external O-rings, where applicable	✓
Leak test of all channels and full functional inspection	✓
Calibration with certificate	✓
Maintenance and repair Available for INTEGRA pipettes with a serial number higher than 7000000	If required
Option: As-found validation of performance before servicing, offering traceability As-found calibration with certificate	Available on request

Please note:

- For on-site calibration, at least 30 pipettes are required per service call.
- For additional info about available calibration certification or measurement methods, please contact your local INTEGRA representative.



INTEGRA Biosciences AG
7205 Zizers, Switzerland
T +41 81 286 95 30
F +41 81 286 95 33
info@integra-biosciences.com

INTEGRA Biosciences Corp.
Hudson, NH 03051, USA
T +1 603 578 5800
F +1 603 577 5529
info-us@integra-biosciences.com

INTEGRA Biosciences Ltd.
Thatcham, Berks RG19 4EP, UK
T: +44 1635 797000
F: +44 1635 797001
info-uk@integra-biosciences.com



www.integra-biosciences.com