

## Overview: Supported inkjet printers (1)

| Supported printers:                                                                                                                                                                                                                                                                  | MI 9330                                                                                                                                               | MI 9029                                                                                                                                                | MI 9410                                                                                                                                                 | MI 9450                                                                                                                                                 |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>MI 9029</li> <li>MI 9330</li> <li>MI 9410</li> <li>MI 9450</li> <li>MI 9018 (discont.)</li> <li>MI 9020 (discont.)</li> <li>MI 9028 (discont.)</li> <li>MI 9030 (discont.)</li> <li>MI 9232 (discont.)</li> <li>MI 9341 (discont.)</li> </ul> | <br><a href="http://www.markem-imaje.com/">www.markem-imaje.com/</a> | <br><a href="http://www.markem-imaje.com/">www.markem-imaje.com/</a> | <br><a href="http://www.markem-imaje.com/">www.markem-imaje.com/</a> | <br><a href="http://www.markem-imaje.com/">www.markem-imaje.com/</a> |
| Performance                                                                                                                                                                                                                                                                          | 4 Lines<br>Alphanumeric<br>wide choice of characters<br>Graphics/Logo printing<br>Side wall printing                                                  | 4 Lines<br>Alphanumeric<br>wide choice of characters<br>Graphics/Logo printing<br>Side wall printing                                                   | 5 Lines<br>Alphanumeric, bar codes,<br>2D codes<br>Graphics/Logo printing<br>Side wall printing                                                         | 5 Lines<br>Alphanumeric, bar codes,<br>2D codes<br>Graphics/Logo printing<br>Side wall printing                                                         |
| Other features                                                                                                                                                                                                                                                                       | Bench top<br>Touch screen<br>USB port<br>RS-232 / Ethernet                                                                                            | Bench top<br>Key pad<br>USB port<br>RS-232                                                                                                             | Bench top<br>Touch screen<br>USB and SD port<br>Ethernet                                                                                                | Bench top<br>Touch screen<br>USB and SD port<br>Ethernet                                                                                                |
| Dimensions (HxWxD)                                                                                                                                                                                                                                                                   | 595 x 303 x 280 mm                                                                                                                                    | 595 x 307 x ??? mm                                                                                                                                     | 555 x 360 x 357 mm                                                                                                                                      | 555 x 360 x 357 mm                                                                                                                                      |
| Weight (kg)                                                                                                                                                                                                                                                                          | 18.5                                                                                                                                                  | 18.5                                                                                                                                                   | 25                                                                                                                                                      | 25                                                                                                                                                      |
| List Price (EUR)                                                                                                                                                                                                                                                                     | <10'000                                                                                                                                               | <10'000                                                                                                                                                | >10'000                                                                                                                                                 | >10'000                                                                                                                                                 |
| Maintenance                                                                                                                                                                                                                                                                          | ink exchange in < 6 min<br>ink runs > 8000 hours                                                                                                      | ink exchange in < 6 min<br>ink runs > 8000 hours<br>automatic cleaning system                                                                          | Maintenance: every 18 month<br>Head cleaning recommended<br>every second week<br>Head can be left attached on<br>MEDIAJET                               | Maintenance: every 18 month<br>Head cleaning recommended<br>every second week<br>Head can be left attached on<br>MEDIAJET                               |
| INTEGRA comments                                                                                                                                                                                                                                                                     | <b>Recommended MI printer,</b><br>lower consumable and<br>maintenance costs.<br>Industrial interface must be<br>ordered.                              | Industrial interface must be<br>ordered.                                                                                                               | Industrial interface must be<br>ordered.                                                                                                                | Industrial interface must be<br>ordered. Higher line speed and<br>storage capacity.                                                                     |

## Overview: Supported inkjet printers (2)

| Supported printers:                                                                                                                                                                                                                        | Linx 89xx<br>(89xx: -00, -10, -20, -40)                                                                                                         | Linx 10                                                                                                                                          | Domino Ax150i                                                                                                                                               | Domino Ax350i                                                                                                                                               |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>Linx 4900, Linx 5900 (printhead MK7)</li> <li>Linx CJ400, Linx 89xx (printhead MK11)</li> <li>Linx 10 (printhead MK13)</li> <li>Domino A320i, Domino A420i, Domino Ax150i, Domino Ax350i</li> </ul> | <br><a href="http://www.linxglobal.com">www.linxglobal.com</a> | <br><a href="http://www.linxglobal.com">www.linxglobal.com</a> | <br><a href="http://www.domino-printing.com">www.domino-printing.com</a> | <br><a href="http://www.domino-printing.com">www.domino-printing.com</a> |
| Performance                                                                                                                                                                                                                                | 6 Lines<br>Alphanumeric, barcodes, 2D codes<br>Graphics/Logo 25 px<br>Side wall / bottom printing                                               | 3 lines<br>Code height 1.8 – 8.8 mm<br>Upper/lower case, graphics<br>Side wall / bottom printing                                                 | 3 Lines<br>Alphanumeric, bar codes, 2D codes<br>Graphics/Logo printing<br>Side wall printing                                                                | 5 Lines<br>Alphanumeric, bar codes, 2D codes<br>Graphics/Logo printing<br>Side wall printing                                                                |
| Other features                                                                                                                                                                                                                             | Benchtop<br>Touch screen<br>USB port<br>Ethernet                                                                                                | Benchtop<br>Touch screen<br>Remote connectivity option<br>Quick and easy to set up                                                               | Bench top<br>Touch screen<br>Optional: USB port, Ethernet, RS232, GPIO                                                                                      | Bench top<br>Touch screen<br>Optional: USB port, Ethernet, RS232, GPIO                                                                                      |
| Dimensions (HxWxD)                                                                                                                                                                                                                         | 503 x 440 x 355 mm                                                                                                                              | 228 x 551 x 189 mm                                                                                                                               | 373 x 400 x 393 mm                                                                                                                                          | 411 x 430 x 331 mm                                                                                                                                          |
| Weight (kg)                                                                                                                                                                                                                                | 24                                                                                                                                              | 11                                                                                                                                               | 16                                                                                                                                                          | 22                                                                                                                                                          |
| List Price (EUR)                                                                                                                                                                                                                           | >10'000                                                                                                                                         | < 8'000                                                                                                                                          | 8'300                                                                                                                                                       | >10'000                                                                                                                                                     |
| Maintenance                                                                                                                                                                                                                                | Maintenance: every 24 months<br>Head cleaning every 3 months                                                                                    | Service module (Easi-Change) to be exchanged after 18 months or 4.000 h / self-service possible                                                  | after 3500 working hours no external service necessary, 10 min. maintenance<br>Head can be left attached on MEDIAJET                                        | after 3500 working hours no external service necessary, 10 min. maintenance<br>Head can be left attached on MEDIAJET                                        |
| INTEGRA comments                                                                                                                                                                                                                           | Recommended, if 2D barcodes are needed. Reliable.                                                                                               | <b>Recommended Linx printer</b><br>Portable, lightweight, small footprint.                                                                       | <b>Recommended Domino printer</b><br>Recommended, if bar codes or 2D codes needed.                                                                          | Recommended, if bar codes or 2D codes needed.                                                                                                               |

## Ordering information MEDIAJET inkjet printer connections kits

| Supported printers: <ul style="list-style-type: none"><li>MI 9018, MI 9020, MI 9028, MI 9029, MI 9030, MI9330</li><li>MI 9232, MI 9410, MI 9450</li><li>Linx 4900, Linx 5900, Linx 89xx (print head MK7)</li><li>Linx CJ400 (print head MK11)</li><li>Linx 10 (print head MK13)</li><li>Domino A320i, Domino A420i, Domino Ax150i, Domino Ax350i</li></ul> | Markem-Imaje (MI)<br>9018-9030, 9330                                                                                                                                              | Markem-Imaje (MI)<br>9232/94xx                                                      | Linx<br>side wall printing<br>(with print head<br>MK7, MK11, MK13)                  | Linx<br>bottom printing<br>(with print head MK7,<br>MK11, MK13)                     | Domino                                                                              |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|                                                                                                                                                                                                                                                                                                                                                            |                                                                                                  |  |  |  |  |
|                                                                                                                                                                                                                                                                                                                                                            | Part number                                                                                                                                                                       |                                                                                     |                                                                                     |                                                                                     |                                                                                     |
| <b>Printer connection kit:</b><br>with printer head<br>fixing device and<br>interface cable                                                                                                                                                                               | 103080                                                                                                                                                                            | 103091                                                                              | 113845                                                                              | 113840                                                                              | 113060                                                                              |
| <b>Operating instructions</b>                                                                                                                                                                                                                                                                                                                              | Operating instructions for connection kits are included and also available in the download center on <a href="http://www.integra-biosciences.com">www.integra-biosciences.com</a> |                                                                                     |                                                                                     |                                                                                     |                                                                                     |
|                                                                                                                                                                                                                                                                                                                                                            | 103983                                                                                                                                                                            | 103980                                                                              | 103956                                                                              | 103981                                                                              | 103982                                                                              |
| <b>Cart</b><br>(W x D x H (cm):<br>70 x 50 x 65-90)                                                                                                                                                                                                                     | -                                                                                                                                                                                 | -                                                                                   | -                                                                                   | 113841                                                                              | -                                                                                   |
| <b>INTEGRA comments</b>                                                                                                                                                                                                                                                                                                                                    | Compatible with:<br>Ø 90 mm, Ø 60 mm or<br>Ø 35 mm Petri dishes.                                                                                                                  | Compatible with:<br>Ø 90 mm or Ø 60 mm<br>Petri dishes.                             | Compatible with:<br>Ø 90 mm or Ø 60 mm<br>Petri dishes.                             | Compatible with:<br>Ø 90 mm, Ø 60 mm or<br>Ø 35 mm Petri dishes.                    | Compatible with:<br>Ø 90 mm or Ø 60 mm<br>Petri dishes.                             |
| <b>Compatibility with MEDIAJET<br/>(manufactured)</b>                                                                                                                                                                                                                                                                                                      | After Dec 2005 / Feb<br>2007; (SN above:<br>MJ : 10305000 /<br>MJ vario : 11307000)                                                                                               | After Feb 2013;<br>(SN above:<br>MJ : 10313568 /<br>MJ vario : 11313247)            | After Feb 2013;<br>(SN above:<br>MJ : 10313568 /<br>MJ vario : 11313247)            | After Mai 2010;<br>(SN above:<br>MJ : 10310371 /<br>MJ vario : 11310128)            | After Feb 2013;<br>(SN above:<br>MJ : 10313568 /<br>MJ vario : 11313247)            |