

SureFlo™ Reagent Reservoirs

The LOWEST reagent waste on the market!

INTEGRA




What is SureFlo™?

Our unique SureFlo™ array is molded into the bottom of each reservoir. This, combined with an inert, hydrophilic surface treatment prevents liquid pooling and spreads reagents evenly across the base of the reservoir, allowing tips to touch the bottom without sealing off or restricting liquid flow. This innovation reduces dead volumes, saving you money on reagents.

Avoid contamination:
Add only the reagent you need to avoid having to pour reagent back into original container



Easy to read printed graduations



Prevent contamination

Inserts double as lids



Reduce plastic waste



Reusable bases



Disposable inserts

Part No.	Product Description	1 Unit	2-4 Units	5+ Units
4371	10ml Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, Polystyrene	89 CAD	80.19 CAD	75.74 CAD
4372	10ml Reservoir, Sterile, SureFlo™, Four Sleeves of 50 per Case, Polystyrene	186 CAD	167.67 CAD	158.36 CAD
4381	25ml Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, Polystyrene	95 CAD	85.05 CAD	80.33 CAD
4382	25ml Reservoir, Sterile, SureFlo™, Four Sleeves of 50 per Case, Polystyrene	196 CAD	176.18 CAD	166.39 CAD
4351	25ml Divided Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, Polystyrene	139 CAD	125.15 CAD	118.19 CAD
4352	25ml Divided Reservoir, Sterile, SureFlo™, Bulk 200 Pack, Polystyrene	289 CAD	260.01 CAD	245.57 CAD
4356	25ml Divided Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, Polypropylene	178 CAD	160.38 CAD	151.47 CAD
4357	25ml Divided Reservoir, Sterile, SureFlo™, Bulk 200 Pack, Polypropylene	371 CAD	334.13 CAD	315.56 CAD
4361	12 Well Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, 3 ml per well, Polystyrene	138 CAD	124 CAD	117 CAD
4366	12 Well Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, 3 ml per well, Polypropylene	176 CAD	159 CAD	150 CAD
4391	100ml Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, Polystyrene	100 CAD	89.91 CAD	84.92 CAD
4392	100ml Reservoir, Sterile, SureFlo™, Four Sleeves of 50 per Case, Polystyrene	203 CAD	182.25 CAD	172.13 CAD
4306	10ml Reservoir Base, 10 Pack	35 CAD	31.59 CAD	29.84 CAD
4304	25ml Reservoir Base, 10 Pack	35 CAD	31.59 CAD	29.84 CAD
4305	100ml Reservoir Base, 10 Pack	35 CAD	31.59 CAD	29.84 CAD

SureFlo™ Reagent Reservoirs


Now In Low Cost, Space Saving Sterile Sleeves




- You asked and we responded! Not everyone needs the large quantity of inserts provided with INTEGRA's bulk pack of reagent reservoirs. Ordering in bulk can be expensive and take up precious lab space.
- INTEGRA is now offering individual sleeves containing only 50 reagent reservoirs. Now every lab can afford to have this product on the lab bench ready to use!

For a limited time only, purchase any of the space saving sleeves listed below and receive a free reservoir base!

Three base versions are available:



10ml



25ml, 25ml-Divided, and 12 Well



100ml

Part No.	Product Description	1 Unit	2-4 Units	5+ Units
4353	Divided Reservoir, 25 ml, Sterile, SureFlo™, 1 Sleeve of 50, Polystyrene	95 CAD	85.05 CAD	80.33 CAD
4358	Divided Reservoir, 25 ml, Sterile, SureFlo™, 1 Sleeve of 50, Polypropylene	120 CAD	108.14 CAD	102.13 CAD
4373	10 ml Reservoir, Sterile, SureFlo™, 1 Sleeve of 50, Polystyrene	61 CAD	54.68 CAD	51.64 CAD
4383	25 ml Reservoir, Sterile, SureFlo™, 1 Sleeve of 50, Polystyrene	63 CAD	57.11 CAD	53.93 CAD
4393	100 ml Reservoir, Sterile, SureFlo™, 1 Sleeve of 50, Polystyrene	66 CAD	59.54 CAD	56.23 CAD



To benefit from this offer, please contact INTEGRA Biosciences by email or phone and mention the promotion code: **SureFloSleeve2021**, or clearly write the promotion code on all faxed purchase orders. This offer is valid until **May 31st, 2021**

INTEGRA Biosciences, US
22 Friars Drive
Hudson, NH 03051

Phone: 1-603-578-5800
Fax: 1-603-577-5529
E-Mail: us@integra-biosciences.com

www.integra-biosciences.com

INTEGRA

Offer valid in Canada only.