

SureFlo™ Reagent Reservoirs

The LOWEST reagent waste on the market!

INTEGRA




What is SureFlo™?

Our unique SureFlo™ array is molded into the bottom of each reservoir. This, combined with an inert, hydrophilic surface treatment prevents liquid pooling and spreads reagents evenly across the base of the reservoir, allowing tips to touch the bottom without sealing off or restricting liquid flow. This innovation reduces dead volumes, saving you money on reagents.

Avoid contamination:
Add only the reagent you need to avoid having to pour reagent back into original container



Easy to read printed graduations



Prevent contamination

Inserts double as lids

Reduce plastic waste



Reusable bases



Disposable inserts

Part No.	Product Description	1 Unit	2-4 Units	5+ Units
4371	10ml Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, Polystyrene	70 CHF	63 CHF	59.50 CHF
4372	10ml Reservoir, Sterile, SureFlo™, Four Sleeves of 50 per Case, Polystyrene	145 CHF	130.50 CHF	123.25 CHF
4381	25ml Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, Polystyrene	70 CHF	63 CHF	59.50 CHF
4382	25ml Reservoir, Sterile, SureFlo™, Four Sleeves of 50 per Case, Polystyrene	145 CHF	130.50 CHF	123.25 CHF
4351	25ml Divided Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, Polystyrene	105 CHF	94.50 CHF	89.25 CHF
4352	25ml Divided Reservoir, Sterile, SureFlo™, Bulk 200 Pack, Polystyrene	218 CHF	196.20 CHF	185.30 CHF
4356	25ml Divided Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, Polypropylene	135 CHF	121.50 CHF	114.75 CHF
4357	25ml Divided Reservoir, Sterile, SureFlo™, Bulk 200 Pack, Polypropylene	280 CHF	252 CHF	238 CHF
4361	12 Well Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, 3 ml per well, Polystyrene	105 CHF	94.50 CHF	89.25 CHF
4366	12 Well Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, 3 ml per well, Polypropylene	135 CHF	121.50 CHF	114.75 CHF
4391	100ml Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, Polystyrene	67 CHF	60.30 CHF	56.95 CHF
4392	100ml Reservoir, Sterile, SureFlo™, Four Sleeves of 50 per Case, Polystyrene	155 CHF	139.50 CHF	131.75 CHF
4306	10ml Reservoir Base, 10 Pack	32 CHF	28.80 CHF	27.20 CHF
4304	25ml Reservoir Base, 10 Pack	32 CHF	28.80 CHF	27.20 CHF
4305	100ml Reservoir Base, 10 Pack	32 CHF	28.80 CHF	27.20 CHF

SureFlo™ Reagent Reservoirs


Now In Low Cost, Space Saving Sterile Sleeves




- You asked and we responded! Not everyone needs the large quantity of inserts provided with INTEGRA's bulk pack of reagent reservoirs. Ordering in bulk can be expensive and take up precious lab space.
- INTEGRA is now offering individual sleeves containing only 50 reagent reservoirs. Now every lab can afford to have this product on the lab bench ready to use!

For a limited time only, purchase any of the space saving sleeves listed below and receive a free reservoir base!

Three base versions are available:



10ml



25ml, 25ml-Divided, and 12 Well



100ml

Part No.	Product Description	1 Unit	2-4 Units	5+ Units
4353	Divided Reservoir, 25 ml, Sterile, SureFlo™, 1 Sleeve of 50, Polystyrene	71 CHF	63.90 CHF	60.35 CHF
4358	Divided Reservoir, 25 ml, Sterile, SureFlo™, 1 Sleeve of 50, Polypropylene	91 CHF	81.90 CHF	77.35 CHF
4373	10 ml Reservoir, Sterile, SureFlo™, 1 Sleeve of 50, Polystyrene	47 CHF	42.30 CHF	39.95 CHF
4383	25 ml Reservoir, Sterile, SureFlo™, 1 Sleeve of 50, Polystyrene	47 CHF	42.30 CHF	39.95 CHF
4393	100 ml Reservoir, Sterile, SureFlo™, 1 Sleeve of 50, Polystyrene	50 CHF	45 CHF	42.50 CHF



To benefit from this offer, please contact INTEGRA Biosciences by email or phone and mention the promotion code: **SureFloSleeve2021**, or clearly write the promotion code on all faxed purchase orders. This offer is valid until **May 31st, 2021**

INTEGRA Biosciences AG
7205 Zizers, Switzerland

T +41 81 286 95 55
F +41 81 286 95 07
info-ch@integra-biosciences.com

www.integra-biosciences.com

INTEGRA

Offer valid in Switzerland only.