

SureFlo™ Reagent Reservoirs

The LOWEST reagent waste on the market!

INTEGRA




What is SureFlo™?

Our unique SureFlo™ array is molded into the bottom of each reservoir. This, combined with an inert, hydrophilic surface treatment prevents liquid pooling and spreads reagents evenly across the base of the reservoir, allowing tips to touch the bottom without sealing off or restricting liquid flow. This innovation reduces dead volumes, saving you money on reagents.

Avoid contamination:
Add only the reagent you need to avoid having to pour reagent back into original container



Easy to read printed graduations



Reduce plastic waste



Reusable bases



Disposable inserts

Prevent contamination

Inserts double as lids



Part No.	Product Description	1 Unit	2-4 Units	5+ Units
4371	10ml Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, Polystyrene	\$66	\$59.40	\$56.10
4372	10ml Reservoir, Sterile, SureFlo™, Four Sleeves of 50 per Case, Polystyrene	\$138	\$124.20	\$117.30
4381	25ml Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, Polystyrene	\$70	\$63.00	\$59.50
4382	25ml Reservoir, Sterile, SureFlo™, Four Sleeves of 50 per Case, Polystyrene	\$145	\$130.50	\$123.25
4351	25ml Divided Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, Polystyrene	\$103	\$92.70	\$87.55
4352	25ml Divided Reservoir, Sterile, SureFlo™, Bulk 200 Pack, Polystyrene	\$214	\$192.60	\$181.90
4356	25ml Divided Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, Polypropylene	\$132	\$118.80	\$112.20
4357	25ml Divided Reservoir, Sterile, SureFlo™, Bulk 200 Pack, Polypropylene	\$275	\$247.50	\$233.75
4361	12 Well Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, 3 ml per well, Polystyrene	\$103	\$92.70	\$87.55
4366	12 Well Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, 3 ml per well, Polypropylene	\$132	\$118.80	\$112.20
4391	100ml Reservoir, Sterile, SureFlo™, 30 Individually Wrapped, Polystyrene	\$74	\$66.60	\$62.90
4392	100ml Reservoir, Sterile, SureFlo™, Four Sleeves of 50 per Case, Polystyrene	\$150	\$135.00	\$127.50
4306	10ml Reservoir Base, 10 Pack	\$26	\$23.40	\$22.10
4304	25ml Reservoir Base, 10 Pack	\$26	\$23.40	\$22.10
4305	100ml Reservoir Base, 10 Pack	\$26	\$23.40	\$22.10

SureFlo™ Reagent Reservoirs


Now In Low Cost, Space Saving Sterile Sleeves




- You asked and we responded! Not everyone needs the large quantity of inserts provided with INTEGRA's bulk pack of reagent reservoirs. Ordering in bulk can be expensive and take up precious lab space.
- INTEGRA is now offering individual sleeves containing only 50 reagent reservoirs. Now every lab can afford to have this product on the lab bench ready to use!

For a limited time only, purchase any of the space saving sleeves listed below and receive a free reservoir base!

Three base versions are available:



10ml



25ml, 25ml-Divided, and 12 Well



100ml

Part No.	Product Description	1 Unit	2-4 Units	5+ Units
4353	Divided Reservoir, 25 ml, Sterile, SureFlo™, 1 Sleeve of 50, Polystyrene	\$70	\$63.00	\$59.50
4358	Divided Reservoir, 25 ml, Sterile, SureFlo™, 1 Sleeve of 50, Polypropylene	\$89	\$80.10	\$75.65
4373	10 ml Reservoir, Sterile, SureFlo™, 1 Sleeve of 50, Polystyrene	\$45	\$40.50	\$38.25
4383	25 ml Reservoir, Sterile, SureFlo™, 1 Sleeve of 50, Polystyrene	\$47	\$42.30	\$39.95
4393	100 ml Reservoir, Sterile, SureFlo™, 1 Sleeve of 50, Polystyrene	\$49	\$44.10	\$41.65



To benefit from this offer, please contact INTEGRA Biosciences by email or phone and mention the promotion code: **SureFloSleeve2021**, or clearly write the promotion code on all faxed purchase orders. This offer is valid until **May 31st, 2021**

INTEGRA Biosciences, US
22 Friars Drive
Hudson, NH 03051

Phone: 1-603-578-5800
Fax: 1-603-577-5529
E-Mail: us@integra-biosciences.com

www.integra-biosciences.com

INTEGRA

Offer valid in United States only.