

	WELLJET dispenser	WELLJET dispenser Stacker
Temperature range	5 - 35°C	
Humidity range	Max. rel. humidity 80% for temperatures up to 31°C, decreasing linearly to 50% rel. humidity at 40°C.	
Dimensions (W x D x H)	20 x 46 x 29 cm (7.9 x 11.4 x 18.1 in)	46 x 46 x 63/102* cm (18.1 x 18.1 x 24.8/40.2* in) * tower for 25/50 plates
Weight	8.8 kg (19.4 lbs)	22.1 kg (48.7 lbs) with 2 towers for 25 plates
Electricity supply	Mains adapter input: 100–240 VAC, 47–63 Hz Device input: 22.8–25.2 VDC, 100 W	
Noise emission	<60 dBA	
User interface	17.8 cm / 7" high resolution touch screen	
Integration	Ethernet interface API commands available	-
Detection method	-	optical sensor
Compatible plate formats	6, 12, 24, 48, 96, 384, 1536, shallow & deep well plates	
Compatible plate heights	5 - 64 mm	6 - 23.5 mm (stacker mode, with lid)
		6 - 45 mm (stacker mode, without lid)
		5 - 64 mm (manual mode)