

PRESS RELEASE

**Pipetting Comfort & Performance Fit For A King ..**

**INTEGRA Biosciences** ([www.integra-biosciences.com](http://www.integra-biosciences.com)) has announced a new addition to their popular PIPETBOY comfort pipetter range – the **PIPETBOY comfort royal**.

The new transparent blue pipetter has been ergonomically designed to fit perfectly into even small hands. Weighing a lightweight 190 grams and powered by 2 NiMH rechargeable batteries the PIPETBOY comfort royal offers scientists a pipetter they can use continuously all day and still remain fatigue-free.

Designed for comfortable single handed pipetting in the 1-50 millilitre volume range the PIPETBOY comfort royal ensures precise and safe operation for applications ranging from rapid pipetting to sensitive precision dosing. A combined all-in-one stand and charging station allows the PIPETBOY comfort royal to be wall mounted or sit on the laboratory bench such that it is always ready to hand. For sensitive samples all parts of the pipetter that come into contact with liquid may be repeatedly autoclaved clean.

Unlike other lower cost pipetting devices the PIPETBOY comfort royal has in-built features to prevent accidental misuse. Incorporating a sterile PTFE-membrane filter the PIPETBOY comfort royal prevents contamination from overflowing as the penetration of liquid is prevented. In addition the PIPETBOY comfort royal will not squirt out under pressure even if a filled pipette is unintentionally aspirated. Incorporating a silicone pipette mount the new pipetter guarantees a sure grip with all commercially available glass and plastic pipettes from 1ml to 50 ml.

For further information on the PIPETBOY comfort royal contact INTEGRA Biosciences on telephone +41-81-286-9530 or email [info@integra-biosciences.com](mailto:info@integra-biosciences.com)

OCTOBER 2004

ibspr43.doc

ILLUSTRATIVE IMAGE:

